

ARDO - Obstetric Products

ARDO - Obstetric Products

Professional - Caring - Proven



The AmniHook® Amniotic Membrane Perforator

The AmniHook® Amniotic Membrane Perforator provides physicians and midwives with a safe and convenient method of rupturing the amniotic membrane at the first attempt. Each AmniHook® is subjected to a rigorous laser-vision inspection, to ensure that only one instrument will be needed per procedure.

- Rounded tip and protected point designed to rupture the amniotic membrane without injury to mother or foetus
- Disposable, sterile and individually packaged for convenience
- 10" long (26.5cm), with moulded sides and a non-slip grip for more precise control

Stock No.	Description	Unit of Sale
65.00.34	sterile Membrane Perforator	100/box



Double-Grip® Umbilical Cord Clamp

An innovative laser-vision inspection system ensures that this precision-engineered cord clamp lives up to its well-earned reputation for excellence. Only clamps which meet our inspection standards 100% will be found in the delivery room.

- Safe, disposable, easy-to-use
- Slipping during use is prevented by means of a non-slip surface on both sides of the clamp
- Maintains a tight seal and exerts constant pressure on any sized cord, even as the cord shrinks
- Hinge-guard confines the cord between serrated jaws to eliminate movement
- No rough edges to snag the clothing of a newborn, or to get caught on bed linen

Stock No.	Description	Unit of Sale
65.00.31	sterile Cord Clamp	125/box
65.00.29	non-sterile Cord Clamp	250/box



Cord Clamp® Clipper

Developed to complement our Double-Grip® Umbilical Cord Clamp, our Cord Clamp® Clipper is specially designed to safely cut through the clamp hinge to remove the clamp from baby.

- Concealed cutting surface eliminates risk of sharp or pointed objects near baby
- Single patient use only, to prevent risk of cross-contamination

Stock No.	Description	Unit of Sale
65.00.32	Single Clipper	1/box



U-Bag® Urine Specimen Collector

The U-Bag® Urine Specimen Collectors are the first choice for convenience, protection, and care for an infant's tender skin. This complete 'system' of urine-collectors offers features which facilitate every infant urine-collection need.

- See-through butterfly patch makes application easy and accurate
- Comfortable for baby, and suitable for use with both males and females
- Construction with a no-flow back valve keeps urine away from the skin to prevent irritation
- Adhesive patch can be folded in half to seal the specimen for sending to the laboratory

U-Bag® with Standard Adhesive Single Specimen Style **sterile**:

Stock No.	Description	Unit of Sale
65.00.07	Newborn (100ml)	100/box
65.00.05	Pediatric (200ml)	100/box

Single Specimen Style **non-sterile**:

Stock No.	Description	Unit of Sale
65.00.04	Newborn (100ml)	100/box
65.00.02	Pediatric (200ml)	100/box



U-Bag® Pt with Porous Cloth Adhesive Single Specimen Style **sterile**:

Stock No.	Description	Unit of Sale
65.00.14	Preemie (10ml)	100/box
65.00.12	Newborn (100ml)	100/box
65.00.10	Pediatric (200ml)	100/box

Single Specimen Style **non-sterile**:

Stock No.	Description	Unit of Sale
65.00.09	Pediatric (200ml)	100/box



U-Bag® Pt with Standard Adhesive 24-hour Style **sterile**:

Stock No.	Description	Unit of Sale
65.00.08	Newborn (100ml)	50/box
65.00.06	Pediatric (200ml)	50/box

24-hour Style **non-sterile**:

Stock No.	Description	Unit of Sale
65.00.03	Pediatric (200ml)	50/box

U-Bag® Pt with Porous Cloth Adhesive 24-hour Style **sterile**:

Stock No.	Description	Unit of Sale
65.00.13	Newborn (100ml)	50/box
65.00.11	Pediatric (200ml)	50/box



PlastiBell® Circumcision Device

A circumcision procedure which requires no special post-operative care or dressing means a significant saving in both time and money. The disposable PlastiBell® Circumcision Device eliminates the need for the repetitive and costly sterilization required for stainless steel clamps. No chance of lost or mis-matched parts that often render other clamps unsafe or unusable.

- Clear, plastic construction provides visual access at all times
- A clean line of excision which helps to promote rapid healing
- No dressings required
- Unique shaped handles for easy size identification
- Sterile, individually packaged

Hospital Assortment Pack:

Stock No.	Description	Unit of Sale
65.00.15	1.1cm (25each), 1.2cm (35each), 1.3cm (45each), 1.4cm (10each), 1.5cm (5each), 1.7cm (5each)	125/box

Physician Office Assortment Pack:

Stock No.	Description	Unit of Sale
65.00.22	1.1cm (5each), 1.2cm (7each), 1.3cm (8each), 1.4cm (3each), 1.5cm (1each), 1.7cm (1each)	25/box

Physician Office Single Sizes (Bulk):

Stock No.	Description	Unit of Sale
65.00.23	1.1cm	25/box
65.00.24	1.2cm	25/box
65.00.25	1.3cm	25/box
65.00.26	1.4cm	25/box
65.00.27	1.5cm	25/box
65.00.28	1.7cm	25/box



Company stamp



Ardo medical Ltd.
Unit 1, Belvedere Trading Estate
Taunton TA1 1BH
United Kingdom

T +44 (0)1823 336 362
F +44 (0)1823 336 364
info@ardomedical.co.uk
www.ardomedical.co.uk

Ardo medical GmbH
Argelsrieder Feld 10
82234 Oberpfaffenhofen
Germany

T +49 (0)8153 90 877 0
F +49 (0)8153 90 877 66
info@ardomedical.de
www.ardomedical.de

Ardo medical AG
Gewerbestrasse 19
6314 Unterägeri
Switzerland

T +41 (0)41 754 70 70
F +41 (0)41 754 70 71
info@ardo.ch
www.ardo.ch