


ARDO - Obstetric Products

Product Materials

Product	Materials
AmniHook®	High impact polystyrene resin, silicone free, PVC-free, latex free
Double-Grip®	Impact modified nylon 66 molding resin, silicone free, PVC-free, latex free
CordClamp®	Handles are molded of Delrin 900 nc10 natural resin from E.I. DuPont. Virgin resin only. Razor blade stock shall be made of 301 stainless steel 45-50 Rockwell C hardness. Eyelet is nickel plated brass. Silicone free, PVC-free, latex free
U-Bag®	12% EVA/LLDPE Blend, silicone free, PVC-free, latex free
PlastiBell®	Cellulose acetate butyrate molding resin. Ligature-Pure flax, vegetable fiber linen thread, 20/4 ply, left twist, 7 turns per inch, polished finish, silicone free, PVC-free, latex free



Ardo medical Ltd. Unit 1, Belvedere Trading Estate Taunton TA1 1BH United Kingdom
 T +44 (0)1823 336 362 F +44 (0)1823 336 364 info@ardomedical.co.uk www.ardomedical.co.uk

Ardo medical GmbH Argelsrieder Feld 10 82234 Oberpfaffenhofen Germany
 T +49 (0)8153 90 877 0 F +49 (0)8153 90 877 66 info@ardomedical.de www.ardomedical.de

Ardo medical AG Gewerbestrasse 19 6314 Unterägeri Switzerland
 T +41 (0)41 754 70 70 F +41 (0)41 754 70 71 info@ardo.ch www.ardo.ch