

# Aspen Surgical® Product Catalog

Blades and Scalpels, Instrument Care, O.R. Accessories,  
Patient and Staff Safety, and Wound Care



Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**



At its foundation, the practice of medicine is all about trust.

Patients put their trust in their physicians; physicians rely on their staff. And all over the world, medical professionals rely on our products as they care for their patients. That is an enormous amount of trust, and we work tirelessly to earn and keep that trust, every day.

We are Aspen Surgical, part of the Hill-Rom Surgical Solutions, and a leading manufacturer of surgical essentials and medical products for practices and providers all over the world.

Aspen Surgical products equip your O.R.'s from floor to ceiling, helping to create a safer and more efficient surgical environment. We are the company that delivers some of the most familiar brands in medicine: the brand behind the brands you trust.

Learn more about all that Aspen Surgical has to offer in the pages of this catalog, or online at [www.aspensurgical.com](http://www.aspensurgical.com).



**Hill-Rom**

Table of Contents

Blades + Scalpels... 8-11

- Bard-Parker® Conventional Blades & Handles pg... 8-9
- Bard-Parker® Protected Blades & Handles pg... 10
- Bard-Parker® Disposable Scalpels (Conventional & Protected) pg... 11

Instrument Care... 14-23

- Oxboro™ Instrument Tip Protectors pg... 14-17
- ColorTrack® Instrument Identification Tape pg... 18-19
- Oxboro™ Corner Protectors pg... 20
- SafeClean® Instrument and Channel Cleaning Brushes pg... 21-23

O.R. Accessories... 26-31

- Anti-Fog, Surgical Headband, Endo Kittner, & Tip Cleaner pg... 26
- Sterion® and Richard-Allan® Vessel Loops pg... 27
- Sterion® Clamp Covers and Suture Aid Booties pg... 28
- Oxboro™ Pin Covers pg... 29
- Richard-Allan® Snares and Safety Pins pg... 29
- Standard and Specialty Disposable Light Handles pg... 30
- O.R. Camera Cover pg... 31

Patient + Staff Safety... 34-41

- Colby™ Fluid Management Products pg... 34-36
- SurgiSafe® O.R. Cord Cover pg... 35
- ErgoSupport® Anti-Fatigue Mat pg... 35
- Colby™ Floor Suction Devices pg... 36
- ViraSorb® Solidifier pg... 36
- Richard-Allan®, Securline®, WriteSite®, Marking Pens pg... 37-39
- Medication Labels pg... 39
- Syringe and Scalpel Handling Systems pg... 40-41

Wound Care... 44-64

- Richard-Allan® Surgical Needles pg... 44-59
- Saber® Bulb Evacuator pg... 60
- 3C™ Collection System pg... 60
- S-Vac® Bulb Evacuator pg... 60
- S-Vac® Fluted Channel Drains pg... 61
- Perforated Drains pg... 62
- Secure-It™ Tube Holders pg... 63
- Transparent Dressing pg... 63
- Sterion® Skin Closures pg... 64

The following brands are registered trademarks of Aspen Surgical Products, Inc:  
AquaNot®, Bard-Parker®, CinchMat®, ColorTrack®, Dr. Fog®, Elvis®, ErgoSupport®,  
EvacMat®, Oxboro®, Rib-Back®, Richard-Allan®, Saber®, SafeClean®, Securline®, Sterion®,  
Surgidyne®, SurgiSafe®, S-Vac®, ViraSorb®, WaterBoom®, WaterBug Quiet Floor Suction®,  
and WriteSite®.

Betadine® is a registered trademark of Purdue Products L.P.  
ChloraPrep® is a registered trademark of CareFusion 2200, Inc.  
Duraprep® is a registered trademark of 3M Company.  
Sterrad® is a Trademark of Advanced Sterilization Products (ASP).  
Superfrost® is a registered trademark of Erie Scientific Company.  
TRUMPF®, Trulight® and iLED® are trademarks of TRUMPF GmbH + Co. KG, Germany.

The following brands are trademarks of Aspen Surgical Products, Inc:  
Colby™, DuoGuard™, InnerVent™, SafeSwitch™, SafetyLock™, Secure-It™, and  
WaterBug™.

© 2016 Aspen Surgical Products, Inc. All rights reserved.







Blades and Scalpels


Physicians and hospitals have relied on the Bard-Parker® brand for more than 100 years. They choose Bard-Parker blades and scalpels for their exceptional precision and quality. Today, Aspen Surgical produces Bard-Parker surgical blades, handles, and disposable scalpels in a broad range of conventional and protected styles.

All Bard-Parker blades and scalpels are made from a proprietary stainless steel composition and a manufacturing process that result an exceptionally strong and precise surgical instrument. To ensure exceptional product quality and prompt delivery, we proudly produce Bard-Parker blades and scalpels in the United States, where our robust manufacturing capacity ensures a consistent, reliable product supply. You have good reason to depend on Bard-Parker.

Bard-Parker® Conventional Blades

Bard-Parker Conventional Blades are both thoroughbreds and workhorses. These are the blades that physicians rely on day in and day out, in the O.R., the ER, labor and delivery, and throughout the medical facility.

- High-quality blades for exceptional control during incisions
- Bard-Parker's exclusive carbon steel Rib-Back® design available for many blades
- Reusable handles available for every blades size
- Foil pouch packaging with a unique inner liner for maximum protection of both the blade and the end user
- Not made with natural rubber latex and intended for single-use



Carbon Steel Blades with Rib-Back® Design

Sterile, individually packaged

Part No.	Size	Qty
371110	10	50/Box, 150/Case
371111	11	50/Box, 150/Case
371112	12	50/Box, 150/Case
371115	15	50/Box, 150/Case
371120	20	50/Box, 150/Case
371121	21	50/Box, 150/Case
371122	22	50/Box, 150/Case
371123	23	50/Box, 150/Case

Stainless Steel Blades

Sterile, individually packaged

Part No.	Size	Qty
371210	10	50/Box, 150/Case
371211	11	50/Box, 150/Case
371212	12	50/Box, 150/Case
371215	15	50/Box, 150/Case
371220	20	50/Box, 150/Case
371221	21	50/Box, 150/Case
371222	22	50/Box, 150/Case
371223	23	50/Box, 150/Case

Carbon Steel Blades with Rib-Back® Design

Non-Sterile, individually packaged (25 Strips of 6)

Part No.	Size	Qty
371310	10	150/Case
371311	11	150/Case
371312	12	150/Case
371315	15	150/Case
371320	20	150/Case
371321	21	150/Case
371322	22	150/Case
371323	23	150/Case
371324	24	150/Case
371325	25	150/Case
371340*	60	150/Case

Stainless Steel, Special Surgeon's Blades

Sterile, individually packaged

Part No.	Size	Qty
371700	10A Plastic & Eye	50/Box, 150/Case
371712	12B Periodontia	50/Box, 150/Case
371716	15C Plastics	50/Box, 150/Case

SafetyLock™ Carbon Steel Blades with Rib-Back® Design

Sterile, individually packaged

Part No.	Size	Qty
371150	10	50/Box, 150/Case
371151	11	50/Box, 150/Case
371152	12	50/Box, 150/Case
371153	15	50/Box, 150/Case
371154	20	50/Box, 150/Case
371155	21	50/Box, 150/Case
371156	22	50/Box, 150/Case
371157	23	50/Box, 150/Case

Bard-Parker® Conventional Handles

From our slender size 7 to our rugged size 4, there's a Bard-Parker® handle suited to most applications and medical professionals.

- Reusable handles available for any blade size
- Not made with natural rubber latex and intended for single-use



Conventional Handles

Use with Bard-Parker® Blades

Non-Sterile

Part No.	Size	Blade Size	Qty
371030	3	10 – 15	5/Box
371031	3L	10 – 15	5/Box
371032	3LA	10 – 15	5/Box
371040	4	20 – 25	5/Box
371041	4L	20 – 25	5/Box
371050	5	10 – 15	5/Box
371060	6	20 – 25	5/Box
371070	7	10 – 15	5/Box
371080	8	60	5/Box
371090	9	10 – 15	5/Box

\*Autopsy, #8 Handle only



Bard-Parker® Protected Blade System

Sharps injuries can occur at any time: prior to, during, or after a procedure. The Bard-Parker protected blade system helps protect your staff from all of them. We designed the system with a similar weight and feel to our convetional products. It even uses the same stainless steel blades, but adds a locking retractable shield that features simple, one-handed activation.

- Similar weight and feel to our convetional scalpel
- Same stainless steel blade as our convetional scalpel
- Locking retractable shield with easy one-handed activation
- Long handle available with select blade sizes
- Not made with natural rubber latex and intended for single-use

Bard-Parker® Protected Blades

Protected Stainless Steel Blades

Sterile, individually packaged

Use with Bard-Parker® Protected Handles

Part No.	Size	Qty
373910	10	50/Box, 150/Case
373911	11	50/Box, 150/Case
373915	15	50/Box, 150/Case
374305	15C	50/Box, 150/Case
373920	20	50/Box, 150/Case
373921	21	50/Box, 150/Case
373922	22	50/Box, 150/Case
373923	23	50/Box, 150/Case




Bard-Parker® Protected Handles

Protected Handles

Use with Bard-Parker® Protected Surgical Blades

Non-Sterile

Part No.	Size	Protected Blade Size	Qty
374030	3	10, 11, 15, 15C	5/Box
374031	3L	10, 11, 15, 15C	5/Box
374040	4	20, 21, 22, 23	5/Box
374041	4L	20, 21, 22, 23	5/Box



Bard-Parker® SafeSwitch™ Reusable Scalpel Handle

Use with any Surgical Blade

Non-Sterile

Part No.	Size	Qty
ST-1013NS	3	5
ST-1013LNS	3L	5
ST-1014NS	4	5
ST-1017NS	7	5



Bard-Parker® Conventional Disposable Scalpels

Disposable Scalpels deliver the quality and precision of Bard-Parker® in a cost-effective scalpel designed for single-use.

- High quality stainless steel blade
- Available in most common blade sizes
- Not made with natural rubber latex and intended for single-use



Disposable Scalpels Stainless Steel Blade

Sterile, individually packaged

Part No.	Size	Qty
371610	10	10/Box, 100/Case
371611	11	10/Box, 100/Case
371612	12	10/Box, 100/Case
371613	12B	10/Box, 100/Case
371615	15	10/Box, 100/Case
371620	20	10/Box, 100/Case
371621	21	10/Box, 100/Case
371623	23	10/Box, 100/Case

Disposable Scalpels Stainless Steel Blade

Non-Sterile

Part No.	Size	Qty
371630	10	100/Bag, 500/Case
371631	11	100/Bag, 500/Case
371633	15	100/Bag, 500/Case
371634	20	100/Bag, 500/Case
371635	21	100/Bag, 500/Case



Disposable Mini Scalpel Stainless Steel

Blade

Non-Sterile

Part No.	Size	Qty
371641	11	100/Bag, 1000/Case

Bard-Parker® Protected Disposable Scalpels

Like our Protected Blade Sytem, Bard-Parker® Protected Disposable Scalpels provide protection from sharps injuries, but are designed for single-use.

- Similar weight and feel to our conventional disposable scalpel
- Same stainless steel blade as our conventional disposable scalpel
- Locking retractable shield with easy one-handed activation
- Long handle available with select blades sizes
- Not made with natural rubber latex and intended for single-use



Protected Disposable Scalpels Stainless Steel Blade

Sterile, individually packaged

Part No.	Size	Qty
372610	10	10/Box, 100/Case
372611	11	10/Box, 100/Case
372612	12	10/Box, 100/Case
372613	12B	10/Box, 100/Case
372615	15	10/Box, 100/Case
372620	20	10/Box, 100/Case
372621	21	10/Box, 100/Case
372622	22	10/Box, 100/Case
372623	23	10/Box, 100/Case
372710	Size 10, Long Handle	10/Box, 100/Case
372711	Size 11, Long Handle	10/Box, 100/Case
372715	Size 15, Long Handle	10/Box, 100/Case

Protected Disposable Scalpels Stainless Steel Blade

Non-Sterile

Part No.	Size	Qty
373810	10	10 Strip, 50/Case
373811	11	10 Strip, 50/Case
373815	15	10 Strip, 50/Case
373820	20	10 Strip, 50/Case
373821	21	10 Strip, 50/Case
373822	22	10 Strip, 50/Case
373823	23	10 Strip, 50/Case








Oxboro® Instrument Tip Protectors

- Designed to protect valuable instruments from breaking and chipping during sterilization and handling
- Help eliminate reprocessing costs caused by torn pouches or wrap and protect staff from sharp sticks during handling

Osteotome Tip Protectors

- Help Protect Osteotome instrument edges from chipping with these tip protectors
- Available tinted or clear and vented or non-vented to comply with user preferences
- All types allow visual instrument inspection and are gas and steam permeable
- Not made with natural rubber latex and intended for single-use



Clear				Blue			
Part No.	Part No.	Size	Qty/Bag	Part No.	Part No.	Size	Qty/Bag
Vented	Non-Vented	mm		Vented	Non-Vented	mm	50
094031BBG	094021BBG	6.4 x 19	50	094011BBG	094001BBG	6.4 x 19	50
094032BBG	094022BBG	9.5 x 25	50	094012BBG	094002BBG	9.5 x 25	50
094033BBG	094023BBG	12.7 x 25	50	094013BBG	094003BBG	12.7 x 25	50
094034BBG	094024BBG	19 x 25.4	50	094014BBG	094004BBG	19 x 25.4	50
094035BBG	094025BBG	25 x 25	50	094015BBG	094005BBG	25 x 25	50
094036BBG	094026BBG	32 x 32	50	094016BBG	094006BBG	32 x 32	50
094037BBG	094027BBG	38 x 38	50	094017BBG	094007BBG	38 x 38	50
094038BBG	094028BBG	44.5 x 32	50	094018BBG	094008BBG	44.5 x 32	50
094040BBG		Assorted	50	094020BBG		Assorted	50

Standard Tip Protector Sizing Guide

Use this guide to determine the standard tip protector on the following page that will best fit your instrument.



Standard Instrument Tip Protectors

- Available solid, tinted or clear and vented or non-vented
- Gas or steam permeable and can be used in flash or conventional autoclave processes
- Solid colors are radiopaque
- Not made with natural rubber latex and intended for single-use



1.6mm x 19mm	Part No.	Qty/Bag
Solid Vented	091011BBG	100
Solid Non-Vented	091001BBG	100
Tinted Vented	093011BBG	100
Tinted Non-Vented	093001BBG	100
Clear Vented	092011BBG	100
Clear Non-Vented	092001BBG	100



2mm x 19mm	Part No.	Qty/Bag
Solid Vented	091012BBG	100
Solid Non-Vented	091002BBG	100
Tinted Vented	093012BBG	100
Tinted Non-Vented	093002BBG	100
Clear Vented	092012BBG	100
Clear Non-Vented	092002BBG	100



2.8mm x 19mm	Part No.	Qty/Bag
Solid Vented	091013BBG	100
Solid Non-Vented	091003BBG	100
Tinted Vented	093013BBG	100
Tinted Non-Vented	093003BBG	100
Clear Vented	092013BBG	100
Clear Non-Vented	092003BBG	100



3.2mm x 25.4mm	Part No.	Qty/Bag
Solid Vented	091014BBG	100
Solid Non-Vented	091004BBG	100
Tinted Vented	093014BBG	100
Tinted Non-Vented	093004BBG	100
Clear Vented	092014BBG	100
Clear Non-Vented	092004BBG	100



5mm x 25mm	Part No.	Qty/Bag
Solid Vented	091015BBG	100
Solid Non-Vented	091005BBG	100
Tinted Vented	093015BBG	100
Tinted Non-Vented	093005BBG	100
Clear Vented	092015BBG	100
Clear Non-Vented	092005BBG	100



10mm x 19mm	Part No.	Qty/Bag
Solid Vented	091016BBG	100
Solid Non-Vented	091006BBG	100
Tinted Vented	093016BBG	100
Tinted Non-Vented	093006BBG	100
Clear Vented	092016BBG	100
Clear Non-Vented	092006BBG	100



2mm x 9mm x 25mm	Part No.	Qty/Bag
Solid Vented	091017BBG	100
Solid Non-Vented	091007BBG	100
Tinted Vented	093017BBG	100
Tinted Non-Vented	093007BBG	100
Clear Vented	092017BBG	100
Clear Non-Vented	092007BBG	100



2mm x 16mm x 25mm	Part No.	Qty/Bag
Solid Vented	091018BBG	100
Solid Non-Vented	091008BBG	100
Tinted Vented	093018BBG	100
Tinted Non-Vented	093008BBG	100
Clear Vented	092018BBG	100
Clear Non-Vented	092008BBG	100



3mm x 25mm x 25mm	Part No.	Qty/Bag
Solid Vented	091019BBG	100
Solid Non-Vented	091009BBG	100
Tinted Vented	093019BBG	100
Tinted Non-Vented	093009BBG	100
Clear Vented	092019BBG	100
Clear Non-Vented	092009BBG	100

Assorted – Solid	Part No.	Qty/Bag
Solid Vented	091020BBG	100
Solid Non-Vented	091010BBG	100

Assorted – Tinted	Part No.	Qty/Bag
Tinted Vented	093020BBG	100
Tinted Non-Vented	093010BBG	100

Assorted – Clear	Part No.	Qty/Bag
Clear Vented	092020BBG	100
Clear Non-Vented	092010BBG	100



DuoGuard™ Instrument Tip Protectors

- Ensure that your instruments stay in the proper open position during processing with DuoGuard™ Instrument Tip Protectors
- Available vented and non-vented, they are tinted for easy viewing of instrument tips
- Not made with natural rubber latex and intended for single-use




Tinted

Part No. Vented	Part No. Non-Vented	Color	Size mm	Qty
093031bbg	093021BBG	White	6.4 x 19	50
093032bbg	093022BBG	Blue	9.5 x 20	50
093033bbg	093023BBG	Green	12.7 x 20	50
093034bbg	093024BBG	Red	19 x 25.5	50
093035bbg	093025BBG	Yellow	15 x 25	50
093037bbg	093027BBG	Orange	25 x 25	50
093040bbg		Assorted	Assorted	50

Low Temperature Web Guards

- An Aspen exclusive
- These unique guards with a flow through, open weave are designed to be used only with low heat, plasma sterilization processes, such as the Sterrad™ Sterilization System
- Great for helping to protect endoscopes and their accessories
- Not made with natural rubber latex and intended for single-use



Part No.	Color	mm	Qty/Bag
099001BBG	Orange Mesh	10 x 76	24
099002BBG	Green Mesh	15 x 76	24
099003BBG	Blue Mesh	20 x 76	24
099004BBG	Red Mesh	40 x 76	24
099005BBG	Black Mesh	60 x 76	16

InnerVent™ Instrument Tip Protectors

- Offer the performance of vented tip protectors with the safety of non-vented tip protectors
- Raised inner nubs allow improved flow through of steam and gas while preventing the tip protectors from “vacuuming” onto the instrument during drying
- Not made with natural rubber latex and intended for single-use



Tinted



Solid/Radiopaque



Clear



Solid, Clear, Tinted

Part No.	Color	Size mm	Qty/Bag
091117BBG	Orange, Solid	2 x 9 x 25	100
091118BBG	Brown, Solid	2 x 16 x 25	100
091119BBG	Black, Solid	3 x 25 x 25	100
092117BBG	Clear	2 x 9 x 25	100
092118BBG	Clear	2 x 16 x 25	100
092119BBG	Clear	3 x 25 x 25	100
093117BBG	Orange, Tinted	2 x 9 x 25	100
093118BBG	Brown, Tinted	2 x 16 x 25	100
093119BBG	Black, Tinted	3 x 25 x 25	100

Gelpi Retractor Tip Protectors

- Protect expensive Gelpi retractors with Aspen Surgical's exclusive Gelpi Retractor Tip Protectors
- Available in three sizes to fit the most commonly used Gelpi retractors
- Gas and steam sterilization compatible
- Not made with natural rubber latex and intended for single-use



Part No.	Color	Size mm	Qty/Bag
096109BBG	Clear	38 Length	55
096129BBG	Clear	51 Length	45
096139BBG	Clear	76 Length	40

Weitlaner Instrument Tip Protectors

- Help Protect your expensive Weitlaner instruments and reduce reprocessing costs associated with compromised wrap
- Specifically designed for use with up to a 4" Weitlaner retractor and offer easy one-handed application
- Not made with natural rubber latex and intended for single-use



Vented, Solid or Clear

Part No.	Size	Color	Qty/Bag
096409BBG	Small	Magenta	50
096410BBG	Small	Clear	50
096419BBG	Medium	Green	50
096420BBG	Medium	Clear	50

Laparoscopic MIS Tip Protectors

- Venting of these specialty tip protectors provides extra flow through to Trocar probes and Veress needles during high level disinfection or sterilization
- Not made with natural rubber latex and intended for single-use



Vented and Clear

Part No.	Size mm	Applications	Qty/Bag
097011BBG	19	Use for Veress Needle	50
097012BBG	20 x 5	Use for Trocar Probe	25
097013BBG	43 x 5	Use for Trocar Probe	20
097014BBG	20 x 10	Use for Trocar Probe	25
097015BBG	57 x 10	Use for Trocar Probe	20

Endoscope Guards

- Help safeguard the eye pieces of delicate endoscopic equipment during storage with these soft, conformable guards
- Not made with natural rubber latex and intended for single-use

Blue, Non-Vented

Part No.	Description	Qty/Bag
095016BBG	Small Eye Piece	6



Round Tip Protectors

- Designed to help protect valuable instruments, ranging in sizes from the smallest Veress Needle to the largest endoscope, during handling and storage
- Tinted for easy visual inspection and vented for improved steam, gas, flow
- Not made with natural rubber latex and intended for single-use



Vented and Tinted

Part No.	Part No.	Color	Fits Tip Sizes	Qty/Bag
096001BBG	096021BBG	White	2-3	70
096002BBG	096022BBG	Blue	3-4	70
096003BBG	096023BBG	Green	4-5	70
096004BBG	096024BBG	Red	5-6	70
096005BBG	096025BBG	Yellow	6-7	70
096006BBG	096026BBG	Maroon	7-8	70
096007BBG	096027BBG	Orange	8-9	70
096008BBG	096028BBG	Brown	9-10	70
096009BBG	096029BBG	Black	10-11	70
096013BBG	096033BBG	Green	13-14	50
096015BBG	096035BBG	Yellow	15-16	50

ColorTrack® Instrument Identification Tape

- Helps make organizing and tracking valuable instruments easy
- It is flexible, adheres securely to instruments, and will not chip, flake, peel, or crack
- With 75+ colors, symbols, and character options, the color-coding options are nearly endless
- If a change is needed, simply remove the ColorTrack I.D. Tape and replace with your new color choice
- Not made with natural rubber latex and non-sterile

Sheet Tape

- Each 8½" x 11" sheet has 374 total inches of tape and contains:  
44 each of 2 inch strips  
44 each of 1½ inch strips  
220 each of 1 inch strips

Solid Color  
¼" Height / Various widths

Part No.	Description
152001EEA	White
152002EEA	Royal Blue
152003EEA	Kelly Green
152004EEA	Red
152005EEA	Yellow
152006EEA	Purple
152007EEA	Orange
152008EEA	Brown
152009EEA	Black
152010EEA	Gray
152011EEA	Lavender
152012EEA	Beige
152013EEA	Light Blue
152014EEA	Lime Green
152015EEA	Pink
152017EEA	Plum
152018EEA	Mint Green

Dots  
¼" Height/Variou widths

Part No.	Description
158040EEA	Dots, Fuchsia
158041EEA	Dots, Medical Blue
158042EEA	Dots, Emerald Green
158043EEA	Dots, Violet

Stars  
¼" Height / Various widths

Part No.	Description
158030EEA	Stars, Fuchsia
158031EEA	Stars, Medical Blue
158032EEA	Stars, Emerald Green
158033EEA	Stars, Violet

White Stripe  
¼" Height / Various widths

Part No.	Description
154002EEA	White/Royal Blue
154003EEA	White/Kelly Green
154004EEA	White/Red
154005EEA	White/Yellow
154007EEA	White/Orange
154008EEA	White/Brown
154009EEA	White/Black
154011EEA	White/Lavender
154013EEA	White/Light Blue
154014EEA	White/Lime Green
154015EEA	White/Pink
154017EEA	White/Plum

Diagonal Stripes  
¼" Height/Variou widths

Part No.	Description
159001EEA	Diagonal, Purple/Black
159002EEA	Diagonal, Orange/Black
159003EEA	Diagonal, Blue/Black
159004EEA	Diagonal, Rubine/Black
159005EEA	Diagonal, Yellow/Black
159006EEA	Diagonal, Green/Black
159007EEA	Diagonal, White/Black

Tri-Color Stripe  
¼" Height / Various widths

Part No.	Description
155051EEA	Red/White/Blue
155052EEA	Red/White/Black
155053EEA	Red/Yellow/Black
155054EEA	Red/Black/Green

Black Stripe  
¼" Height / Various widths

Part No.	Description
153001EEA	Black/White
153003EEA	Black/Kelly Green
153004EEA	Black/Red
153005EEA	Black/Yellow
153007EEA	Black/Orange
153010EEA	Black/Gray
153011EEA	Black/Lavender
153012EEA	Black/Beige
153013EEA	Black/Light Blue
153014EEA	Black/Lime Green
153015EEA	Black/Pink
153017EEA	Black/Plum

Triangles  
¼" Height/Variou widths

Part No.	Description
158050EEA	Triangles, Fuchsia
158051EEA	Triangles, Medical Blue
158052EEA	Triangles, Emerald Green
158053EEA	Triangles, Violet

Numbers  
¼" Height / Various widths

Each sheet contains numbers 0–8

Part No.	Description
11111	Numbers, Fuchsia
22222	Numbers, Medical Blue
33333	Numbers, Emerald Green
44444	Numbers, Violet

Assorted Symbols  
¼" Height / Various widths

Part No.	Description
158060EEA	Red Hearts on White
158080EEA	Alphabet, Fuschia/Wt, A-E

4-Color  
¼" Height / Various widths

Note: Assorted color package will contain separate sheets each with its own color.

Part No.	Description
152031EEA	Red/Yellow/Purple/Orange
152032EEA	Blue/White/Black/Green

Specialty Stripe  
¼" Height / Various widths

Part No.	Description
156003EEA	Blue/Green
156004EEA	Blue/Red
157002EEA	Yellow/Blue
157003EEA	Yellow/Green
157004EEA	Yellow/Red

Code to Specialty  
¼" Height / Various widths

Part No.	Description
ENT	Mint/White ENT
ER ER	Purple/White ER
GEN	Royal Blue/White GEN
OB/GYN	Lt Blue/White OB/GYN
ORTHO	Kelly Green/White ORTHO
NEURO	Orange/White NEURO
CARDIAC	Red/ White CARDIAC
6	Red/White ASPEN SYMBOL

Roll Tape  
Roll Tape ¼" x 300"

Part No.	Description
151001EEA	White
151002EEA	Royal Blue
151003EEA	Kelly Green
151004EEA	Red
151005EEA	Yellow
151006EEA	Purple
151007EEA	Orange
151008EEA	Brown
151009EEA	Black

Roll Tape ¼" x 250"

Part No.	Description
151011EEA	White
151012EEA	Royal Blue
151013EEA	Kelly Green
151014EEA	Red
151015EEA	Yellow
151016EEA	Purple
151017EEA	Orange
151018EEA	Brown
151019EEA	Black



**Application Instructions:**  
For secure placement through multiple processing cycles, ColorTrack® I.D. Tape should be wrapped 1½ times around the Instrument and excess length should be trimmed. Instrument sterility under ColorTrack® I.D. Tape is assured after processing with properly applied tape.\*

\* Test Results Available On Request  
Size and color of tapes shown may vary from actual products.

Sterion® Same Day Tape Remover

- An easy, instrument-friendly solution to remove ColorTrack® Tape
- Helps eliminate the risk of cuts and sticks common with scalpel or utility knife tape removal
- Neutral scent helps make Same Day Tape Remover comfortable to use

Part No.	Description
151066EEA	1 Gallon Container

Helps to  
Easily  
Remove  
Our Tape!





ColorTrack® Write-On ID Labels

- Ideal for customization of trays, storage areas, cabinets or carts, these color-coded labels are formatted on 8½"x 11" sheets with three different label sizes on each sheet
- For best results use with permanent markers only
- Not made with natural rubber latex



Part No.	Description	Qty/Sheet
161002EEA	Blue Border	44
161003EEA	Green Border	44
161004EEA	Red Border	44
161005EEA	Yellow Border	44

Oxboro® Corner Protectors

- Protect your instrument trays from rips, holes, and tears in sterilization wrap due to sharp tray corners, improper handling, or heavy surgical trays
- Aspen Surgical's Corner Protectors are designed to create a flexible barrier between the tray's sharp corners and the sterilization wrap
- They can be used with steam, gas, ozone, and low temperature plasma sterilization and are reusable
- Not made with natural rubber latex



Part No.	Description	Qty/Bag
CP-80NS	Corner Protectors	80

Sterion® Color-Coded Identification Locking Tags

- These tamper evident locking tags can be customized to a variety of uses such as identifying instruments, carts, trays, containers or cabinets
- Available pre-printed with "Repair", "Sharpen", or "Good" on one side and writable on the other or blank on one side and writable on the other
- Not made with natural rubber latex



Part No.	Description	Qty/Bag
251001BBG	Red, Imprint: Repair, Write-On	100
251002BBG	Yellow, Imprint: Sharpen, Write-On	100
251003BBG	Green, Imprint: Good, Write-On	100
251004BBG	Blue, Write-On	100
251005BBG	White, Write-On	100
251006BBG	Green, Write-On	100
251007BBG	Black, Write-On	100
251008BBG	Purple, Write-On	100
251009BBG	Orange, Write-On	100
251010BBG	Assorted Colors, Write-On	100




SafeClean® Instrument and Channel Cleaning Brushes

- Designed to clean surgical/medical instruments and endoscopic devices
- The complete line of SafeClean® Cleaning Brushes are strong enough to effectively remove surgical residue, yet gentle enough to use on the most delicate instruments and devices
- Not made with natural rubber latex

General Purpose Brushes

- Can be used throughout all of your instrument and device reprocessing areas to clean uniquely shaped and larger bore surgical, medical, and endoscopic instruments and devices
- Not made with natural rubber latex



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/Bag
241007BBG	1mm (0.039")	18"	0.5"	2
241010BBG	1.6mm (0.063")	18"	0.5"	2
241015BBG	1.8mm (0.071")	18"	0.5"	2
241019BBG	2mm (0.079")	12"	1"	2
241020BBG	2.2mm (0.086")	18"	0.5"	2
241024BBG	2.5mm (0.098")	24"	2.5"	2
241025BBG	2.7mm (0.106")	18"	1"	2
241029BBG	3mm (0.118")	12"	1"	2
241030BBG	3.2mm (0.126")	18"	2"	2
241032BBG	4mm (0.158")	12"	1"	2
241033BBG	5mm (0.197")	12"	2.5"	2
241034BBG	5mm (0.197")	24"	2.5"	2
241035BBG	5mm (0.197")	16"	2.5"	2
241040BBG	6mm (0.236")	16"	2.5"	2
241041BBG	6mm (0.236")	24"	2.5"	2
241042BBG	7mm (0.275")	16"	2.5"	2
241043BBG	7mm (0.275")	24"	2.5"	2
241044BBG	10mm (0.393")	16"	2.5"	2
241045BBG	11mm (0.433")	16"	2.5"	2
241046BBG	11mm (0.433")	24"	2.5"	2
241047BBG	10mm (0.393")	24"	2.5"	2
241050BBG	12mm (0.472")	16"	2.5"	2
241055BBG	15mm (0.590")	16"	2.5"	2
241060BBG	20mm (0.787")	16"	2.5"	2

Kits

- Help make ordering easier with all of the brushes you need in one kit

Small Cleaning Brush Kit

- Includes various small general purpose cleaning brushes
- Not made with natural rubber latex

Part No.
241100bbg – 1 Each
241010BBG
241015BBG
241020BBG
241025BBG
241030BBG

Large Cleaning Brush Kit

- Includes various large general purpose cleaning brushes
- Not made with natural rubber latex

Part No.
241200bbg – 1 Each
241035BBG
241040BBG
241045BBG
241050BBG
241055BBG
241060BBG

“Best Sellers” Cleaning Brush Kit

- Includes a selection of our most popular general purpose cleaning brushes
- Not made with natural rubber latex

Part No.
241300bbg – 1 Each
241010BBG
241015BBG
241020BBG
241025BBG
241030BBG
241035BBG
241040BBG
241045BBG
241050BBG
241055BBG
241060BBG

General Cleaning Brush Kit

- Contains an extensive selection of general purpose brushes, as well as toothbrush and fan tip brushes
- Not made with natural rubber latex

Part No.
241600bbg – 3 Each
241025BBG
241030BBG
241035BBG
241040BBG
241041BBG
241045BBG
241046BBG
241050BBG
241055BBG
241060BBG
242001BBG
242002BBG
242003BBG
243001BBG
243004BBG

Micro/Ophthalmic Brush Kit

- Contains a selection of our smallest general purpose cleaning brushes
- Not made with natural rubber latex

Part No.
241700bbg – 3 Each
241007BBG
241010BBG
241015BBG
241020BBG
241025BBG
241030BBG


Tooth Brush Kit

- Not made with natural rubber latex

Part No.
243100bbg – 1 Each
243001BBG
243002BBG
243003BBG

Suction Tube Cleaning Brushes

- Brushes are designed to fit into and effectively clean very small bore devices and instruments used in surgery, MIS, and therapy procedures today
- Not made with natural rubber latex



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
241003BBG	3 French (1mm, 0.039")	6"	0.5"	2
241005BBG	5 French (1.67mm, 0.066")	6"	0.5"	2
241006BBG	6 French (2mm, 0.079")	8"	0.5"	2
241008BBG	7 French (2.3mm, 0.092")	8"	0.5"	2
241009BBG	8 French (2.7mm, 0.105")	8"	0.5"	2
241011BBG	9 French (3mm, 0.118")	8"	0.5"	2
241012BBG	10 French (3.3mm, 0.131")	8"	0.5"	2
241013BBG	11 French (3.7mm, 0.144")	8"	0.5"	2
241014BBG	12 French (4mm, 0.158")	8"	0.5"	2

Fan-Tip Brush

- Features bristles that wrap around the end of the brush to clean bottles, jars, and other large instruments
- Not made with natural rubber latex



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
242001BBG	4.8mm-12.7mm Tapered (0.189"–0.5")	7"	3"	2
242002BBG	13mm (0.512")	6"	3"	2
242003BBG	16mm (0.623")	5"	2.5"	2
242046BBG	9.6mm (0.384")	24"	2.5"	2

Tracheostomy Brushes


- Uniquely designed to clean tubes
- Not made with natural rubber latex



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
247001BBG	Fits Trach Sizes 00–5FR 9mm (0.354")	4.75"	1.85"	2
247002BBG	Fits Trach Sizes 6–12FR 13mm (0.512")	8"	3.66"	2

Toothbrush Style

- Variety of bristle types available
- Designed to clean all types of instruments and stains including nylon, brass, and steel
- Not made with natural rubber latex



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
243001BBG	13mm (0.512"), Nylon	7.25"	1.5"	2
243002BBG	13mm (0.512"), Brass	7.25"	1.5"	2
243003BBG	13mm (0.512"), Steel	7.25"	1.5"	2

Double-Ended Nylon Cleaning Brushes

- Features one end with an extra large brush head for general cleaning and one end with a single row of bristles for more intricate cleaning
- Not made with natural rubber latex




Available Sterile and Non-Sterile

Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty
243004BBG	Non-Sterile			
	Short End, 2mm (0.079")	7"	0.867"	2/Bag
	Long End, 8mm (0.315")		1.38"	
243014PBX	Sterile			
	Short End, 2mm (0.079")	7"	0.867"	1/Pack,
	Long End, 8mm (0.315")		1.38"	12/Packs/Box

Reusable Channel Cleaning Brushes


- Channel Brushes with tapered brush heads, elevator sheath, and comfortable round, plastic handles
- Available in a variety of lengths and materials to fit the range of endoscopic devices
- Not made with natural rubber latex



Part No.	Description	Total Length	Sheath Diameter	Qty/ Bag
248107BBG	3-5mm	70cm	1.6mm	1
248115BBG	3-5mm	115cm	1.6mm	1
248160BBG	3-5mm	160cm	1.6mm	1
248180BBG	3-5mm	230cm	1.8mm	1
248190BBG	3-5mm	120cm	1.8mm	1

Disposable Scope Channel Cleaning Brushes


- Brush heads that range from to 2.5mm - 5mm and will fit scopes as small as 1.1mm
- Features more than double the bristle length than competing low-cost brushes
- Strong, sturdy handle making it easy to pull through the scope, while still cleaning the scope effectively
- Intended for single-use and not made with natural rubber latex



Part No.	Description	Total Length	Sheath Diameter	Qty/ Box
248102BBG	2.5mm	230cm	1mm	50
248103BBG	3mm	230cm	1.55mm	50
248104BBG	5mm – Double Ended	230cm	1.85mm	50
248105BBG	5mm	230cm	1.85mm	50

Hand/Nail Brushes

- Used to clean dirt on hands and nails
- Can also be used for general instrument cleaning
- Not made with natural rubber latex



Part No.	Brush Head Area	Qty/Bag
244001BBG	35mm x 93mm (1.375" x 1.75")	2

Instrument Brush Storage

- Convenient storage with our acrylic racks organize your cleaning brushes
- Helps ensure you will have the right brush for the job
- Not made with natural rubber latex

Part No.	Description	Qty/Box
246001BBG	Brush Rack, Medium (without drip tray)	1
246003BBG	Brush Rack, Large (with drip tray)	1









Dr. Fog® Anti-Fog

- Offered with our easy squeeze bottle and sponge set or as a ready-to-use sponge, which is dry to touch and alcohol free
- Radiopaque
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
DF-3100	Original Solution with Sponge	20
DF-3100B	Original Solution with Sponge	300
DF-3120	Alcohol-Free Dry Treated Sponge	40



Elvis® Anti-Fog

- Offered in both a squeeze bottle and pretreated sponge
- Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Qty/Box
400750	Anti-Fog, Pre-Treated Sponge	10
400751	Anti-Fog Bottle with Sponge	10

Surgical Headband

- Used to wick sweat away from the clinician
- Non-sterile and not made with natural rubber latex




Part No.	Description	Qty/Box
400450	Surgical Headband	10

Endoscopic Kittner


- Our blunt dissecting instrument is 38cm long with a 5mm stainless steel shaft
- Available in both single and dual tip
- Radiopaque tip
- The brushed finish allows for safe use during laser surgery
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
KT-9101	Single Tip	20
KT-9105	Dual Tip	20



Sterion® Cautery Tip Cleaner


- Assists in the cleaning and removing of material from cautery tips during surgical procedures
- Consists of a 1.8" x 1.8" foam pad with abrasive material on one side and strong adhesive on the other for attachment to both disposable and non-disposable drapes
- Cleans both mono-polar and bi-polar probes
- Radiopaque
- Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Qty/Box
032001PBX	Cautery Tip Cleaner	35
TC-100	Cautery Tip Cleaner	100
TC-100A	Cautery Tip Cleaner	25
TC-100B	Cautery Tip Cleaner	240

Sterion® and Richard-Allan® Vessel Loops

- Used in occluding, retracting, and identifying arteries, veins, tendons and nerves in surgical procedures
- Extruded from medical-grade silicone, with the exception of the Richard-Allan® vessel loops, which are made of high performance, medical grade elastomer
- Radiopaque
- Packaged both sterile and non-sterile
- Intended for single-use and not made with natural rubber latex



Sterion® Vessel Loops

2 Loops / Pkg; 10 Pkgs / Box Double Pouched

Part No.	Color	mm Width x Length	Inches Width x Length
011001PBX	Red Mini	0.83 x 406	0.03 x 16
011002PBX	Blue Mini	0.83 x 406	0.03 x 16
011003PBX	Yellow Mini	0.83 x 406	0.03 x 16
011004PBX	White Mini	0.83 x 406	0.03 x 16
011011PBX	Red Maxi	1.0 x 406	0.04 x 16
011012PBX	Blue Maxi	1.0 x 406	0.04 x 16
011013PBX	Yellow Maxi	1.0 x 406	0.04 x 16
011014PBX	White Maxi	1.0 x 406	0.04 x 16
011021PBX	White Supermaxi	1.3 x 559	0.05 x 22
011022PBX	Blue Supermaxi	1.3 x 559	0.05 x 22
011032PBX	Blue Micro	0.41 x 406	0.019 x 16
011301PBX	Red Mini	0.83 x 762	0.03 x 30
011302PBX	Blue Mini	0.83 x 762	0.03 x 30
011303PBX	Yellow Mini	0.83 x 762	0.03 x 30
011304PBX	White Mini	0.83 x 762	0.03 x 30
011311PBX	Red Maxi	1.0 x 762	0.04 x 30
011312PBX	Blue Maxi	1.0 x 762	0.04 x 30
011313PBX	Yellow Maxi	1.0 x 762	0.04 x 30
011314PBX	White Maxi	1.0 x 762	0.04 x 30

Richard-Allan® Vessel Loops



2 Loops / Pkg; 50 Pkgs / Box Banded and Pouched

Part No.	Color	mm Width x Length	Inches Width x Length
3901	Maxi Red	1.0 x 457	0.04 x 18
3902	Maxi Blue	1.0 x 457	0.04 x 18
3903	Maxi White	1.0 x 457	0.04 x 18
3904	Maxi Yellow	1.0 x 457	0.04 x 18
3905	Mini Red	.83 x 457	0.03 x 18
3906	Mini Blue	.83 x 457	0.03 x 18
3907	Mini White	.83 x 457	0.03 x 18
3908	Mini Yellow	.83 x 457	0.03 x 18


4 Loops / Pkg; 25 Pkgs / Box

Part No.	Color	mm Width x Length	Inches Width x Length
3910	Four Color Maxi Set	1.0 x 457	0.04 x 18
3912	Four Color Mini Set	0.83 x 457	0.03 x 18



Sterion® Clamp Covers

- Assist in the atraumatic clamping of delicate vessels during surgery
- Available in radiopaque silicone, non-radiopaque fabric, and radiopaque fabric
- A variety of sizes to provide outstanding convenience and flexibility of use
- Packaged sterile or non-sterile
- Intended for single-use and not made with natural rubber latex




Part No.	Color	Material	Radiopaque	mm Length	Inches Length	Qty/Box
Sterile						
061011PBX	Red	Silicone	Yes	2.5 x 130	0.10 x 5	10 Pkgs of 2
061012PBX	White	Silicone	Yes	5.1 x 130	0.20 x 5	10 Pkgs of 2
061013PBX	Blue	Silicone	Yes	7.6 x 130	0.30 x 5	10 Pkgs of 2
071001PBX	White	Fabric	No	3.2 x 130	0.125 x 5	10 Pkgs of 4
071002PBX	White	Fabric	No	9.5 x 130	0.375 x 5	10 Pkgs of 4
072001PBX	White	Fabric	Yes	3.2 x 150	0.125 x 6	10 Pkgs of 4
072002PBX	White	Fabric	Yes	9.5 x 150	0.375 x 6	10 Pkgs of 4
072004PBX	White	Fabric	Yes	6.3 x 150	0.25 x 6	10 Pkgs of 4

Sterion® Suture Aid Booties

- Protect, tag, and help locate sutures
- Vinyl suture booties are available in two-, three-, or five-pair configurations and are positioned in an adhesive-backed foam block for quick and easy usage
- Radiopaque
- Choose from a variety of sizes and colors as listed below
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Tray	Qty/Box
051003PBX	Standard Yellow	5 Pair	5 Trays
051011PBX	Standard Red	5 Pair	5 Trays
051012PBX	Standard Blue	5 Pair	5 Trays
051014PBX	Standard White	5 Pair	5 Trays
051015PBX	Standard Asst.	5 Pair	5 Trays
051004PBX	Large Pink	5 Pair	5 Trays
053003PBX	Standard Yellow	3 Pair	5 Trays
052004PBX	Large Pink	2 Pair	5 Trays
052003PBX	Standard Yellow	2 Pair	5 Trays
053030PBX	X-tra Mini Red	3 Pair	5 Trays



Oxboro® Pin Covers

- Help to protect patients, staff, and clothing from the sharp ends of Steinmann Pins and K-Wires
- Made of thermoplastic and are available in three color-coded sizes to fit all eleven sizes of Steinmann Pins and K-Wires
- Sterile, intended for single use, and not made with natural rubber latex


Part No.	Color	Fits This Size Range/Inches	Qty/Box
101001PBX	Yellow	0.028 – 0.062	1/ Pouch, 60 Pouches/Box
101002PBX	Blue	0.078 – 0.125	2/Pouch, 30 Pouches/Box
101003PBX	Red	0.141 – 0.188	2/Pouch, 30 Pouches/Box
101000PBX	Assorted	Assorted	50 Assorted Pouches/Box



Richard-Allan® Snares

- Formed in pre-cut, pre-bent sizes, ready for immediate use
- Sterile, intended for single-use, and not made with natural rubber latex


Part No.	Description	Qty/Box
2405	Nasal Snare Size 5	2/Pack, 50 Packs/Box
2406	Nasal Snare Size 6	2/Pack, 50 Packs/Box
2407	Tonsil Snare Size 7	2/Pack, 50 Packs/Box
2408	Tonsil Snare Size 8	2/Pack, 50 Packs/Box
2409	Tonsil Snare Size 9	2/Pack, 50 Packs/Box



Richard-Allan® Safety Pins

- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description Inches	Qty/Box
2201	1	2/Pack, 50 Packs/Box
2202	1.5	2/Pack, 50 Packs/Box
2203	2	2/Pack, 50 Packs/Box



Standard Disposable Light Handles

- Two systems offered for convenience: Rigid or Flexible Light Handle Cover
- Helps eliminate costly reprocessing and create uniformity throughout your hospital
- Light Handle Adaptors are available for most of the popular lights
- Please contact your local Aspen Surgical representative for further adaptor information
- Sterile, intended for single-use, and not made with natural rubber latex


Part No.	Description	Qty/Box
LT-R01	Rigid Handle	24
LT-F01B	Flexible Cover (1 Each)	100
LT-F02	Flexible Cover (2 Each)	80
LT-F03	Flexible Cover (3 Each)	80



Specialty Disposable Light Handle

- Designed to have a similar look and feel as the TRUMPF® reusable handle
- Works with all TRUMPF® ALC plus Surgical lights, including the Trulight® 5000 and iLED® 3/5
- Lens clarity is equivalent to the reusable handle, providing rapid response of the laser and sensor
- Ultra-durable lens and a snap-on feature that helps allow quick and easy installation and removal
- Light intensity can be adjusted with the touch sensors while the disposable handle is being used
- Intended for single-use and not made with natural rubber latex


Part No.	Description	Sterility	Qty/Box
LT-ALC01	Disposable Light Handle	Sterile	18



O.R. Camera Cover

- Helps save time and reduce inconvenience associated with having to sterilize reusable camera covers and helps enable a quicker O.R. turnaround time
- The clear lens allows for complete visualization of the surgical field with no distortion
- Adaptors are currently available for many common O.R. cameras
- Please contact your local Aspen Surgical representative for further adaptor information
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
LT-C02	Camera Cover	20









Colby™ Fluid Management Products


Our Colby™ line provides a convenient and effective means of managing infectious waste fluid generated during wet procedures. Our products help to reduce the likelihood of injuries associated with wet floors. They aid in protecting healthcare professionals and housekeeping personnel from contact with contaminated substances and help support your efforts to comply with OSHA's Bloodborne Pathogens Standard. Our products do not change the contents of medical waste, but rather help your team to more effectively contain it.

SurgiSafe® Absorbent Floor Mats

- Help keep floors clean, dry, and slip-free, even during fluid intense procedures
- Absorb irrigants, blood, and bodily fluids through unique wicking, absorbency, and fluid retention properties, while conforming to uneven floor surfaces without bunching
- Help protect your staff from slips, trips, and falls in the operating room
- Help reduce staff exposure to infectious waste
- Help eliminate re-using blankets and towels on patients
- Help reduce O.R. turnaround time
- Contain fluid in one area for easy clean up
- Non-Sterile, intended for single-use, and not made with natural rubber latex

SurgiSafe® Specialty Absorbent Floor Mats

- Designed to include unique features beyond those offered by the Standard Absorbent Floor Mat Line
- Low profile
- Offered in large sizes
- Blue mat features fluid barrier backing to help keep floor beneath dry and clean
- Yellow mat absorbs fluid from the top and bottom



Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty/Case
82350	Blue	Medium	Yes	23" x 40"	50
84625	Blue	Medium	Yes	46" x 40"	25
84610	Blue	Medium	Yes	46" x 40"	10
84672	Blue	Medium	Yes	40" x 72"	10
81860-2	Yellow	Medium	No	18" x 40"	60
83630-2	Yellow	Medium	No	36" x 40"	30
83610-2	Yellow	Medium	No	36" x 40"	10
83615-2L	Yellow	Medium	No	36" x 72"	15
83625-2L	Yellow	Medium	No	36" x 72"	12
83630-M	Yellow	Medium	No	36" x 40"	15
	White	Medium	No	28" x 40"	15


SurgiSafe® Standard Absorbent Floor Mats

- Made with a minimum of 60% pre-consumer recycled material
- Available in three absorbency levels as well as with and without fluid barrier backing for convenience

**Pink** – High level of fluid absorbency, will absorb 16x's its dry weight

**White** – Medium level of fluid absorbency, will absorb 14x's its dry weight

**Green** – Low level of fluid absorbency, will absorb 11x's its dry weight



Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty/Case
83430	Pink	High	No	28" x 40"	30
83472	Pink	High	No	28" x 72"	30
83330	Pink	High	Yes	28" x 40"	30
83372	Pink	High	Yes	28" x 72"	30
81860	White	Medium	No	18" x 40"	60
818100	White	Medium	No	18" x 18"	100
81850	White	Medium	No	18" x 18"	50
83610	White	Medium	No	28" x 40"	10
83630	White	Medium	No	28" x 40"	30
83625-L	White	Medium	No	28" x 72"	12
83672	White	Medium	No	28" x 72"	30
83530	White	Medium	Yes	28' x 40"	30
83572	White	Medium	Yes	28" x 72"	30
84030	Green	Low	No	28" x 40"	30
84072	Green	Low	No	28" x 72"	30
840125	Green	Low	No	28" x 125"	24
83030	Green	Low	Yes	28" x 40"	30
83072	Green	Low	Yes	28" x 72"	30
830125	Green	Low	Yes	28" x 125"	24

SurgiSafe® CinchMat® Absorbent Floor Mat

- Combines the effective absorbency of the SurgiSafe® mat with an easier way to handle the mat during disposal
- Pull cord design helps make clean-up and O.R. turnaround time simple and fast
- The top layer has wicking abilities that quickly absorbs irrigants, blood, and bodily fluids from the O.R. floor to help prevent pooling and splashing
- Skid-resistant fluid barrier backing helps ensure that the floor beneath the mat remains dry and prevents staining
- Large enough for two people
- Aids in keeping the O.R. floor dry to reduce the risk of injuries related to wet floors
- Non-Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty/Case
84625-C	Blue	Medium	Yes	46" x 40"	10

SurgiSafe® EvacMat® Absorbent Floor Mat

- Combines the functionality of absorption and suction to aid in efficiently removing fluid waste from the O.R. floor
- Extremely quiet operation helps to decrease distraction of the O.R. staff– as quiet as an indoor conversation
- The top layer of the mat contains wicking to help prevent pooling
- Works well with wall suction
- Simple to use and clean up – red bag friendly
- Low profile allows for the O.R. staff to roll equipment over the mat without loss of suction
- Large enough for two people
- Provides cushion to help reduce fatigue during long procedures
- Aids in keeping the O.R. floor dry to reduce the risk of injuries related to wet floors
- Non-Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty/Case
95005	Yellow	Medium	Yes	36" x 48"	5

SurgiSafe® O.R. Cord Cover

- Intended to help organize tubes and cords in the Operating Room
- Serves as a “beacon” so healthcare personnel know where floor cords and tubes are located and helps decrease the chance of staff tripping on them
- Made of a durable material that can be rolled over
- Features a residue-free adhesive, helping it stay in one place
- Non-Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Color	Size	Qty/Case
96050	O.R. Cord Cover	Yellow	9.5" x 24"	50

ErgoSupport® Anti-Fatigue Mat

- Intended to help reduce fatigue that healthcare personnel experience during long hours of standing in one location
- Uses a dual foam and gel technology to provide comfort under the foot
- The top surface provides a necessary friction level with damp or dry surgical booties
- The bottom surface is made of a high-traction material to provide slip resistance
- Available in four sizes
- Water-resistant
- Reusable for up to 6-12 months
- Non-sterile and not made with natural rubber latex



Part No.	Description	Color	Size	Qty/Case
81005	Anti-Fatigue Mat	Red	13" x 15"	5
82001	Anti-Fatigue Mat	Red	18" x 24"	1
82005	Anti-Fatigue Mat	Red	18" x 24"	5
83002	Anti-Fatigue Mat	Red	20" x 48"	2
84003	Anti-Fatigue Mat	Red	20" x 32"	3



Colby™ Floor Suction Devices

WaterBug Quiet Floor Suction® Device


- Efficiently suctions fluid waste to help reduce risks and maintain a dry operating room floor
- Its raised design and 12' of flexible tubing allows for easy repositioning
- Suctions unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Tubing	Qty/Case
90010	WaterBug Quiet Floor Suction® Device	Yes	10
90020	WaterBug Quiet Floor Suction® Device	Yes	20
90030	WaterBug Quiet Floor Suction® Device	Yes	30

AquaNot® Floor Suction Device

- Has similar characteristics as the WaterBug Quiet Floor Suction® Device, but has a lower profile
- Efficiently suctions fluid waste to help reduce risks and maintains a dry operating room floor
- Suctions unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Tubing	Qty/Case
97020	AquaNot® Floor Suction Device	Yes	20
97060	AquaNot® Floor Suction Device	Yes	60
97050-NT	AquaNot® Floor Suction Device	No	50

WaterBoom® Floor Suction Strip

- Creates a dam to help hold back, contain, and quietly suction an unlimited amount of waste fluid
- Provides strong suction, while conforming to your O.R. floor
- Its low profile helps you step over it without suction loss
- Its 12' tubing attaches from the side to help reduce tripping hazards
- This unobtrusive device does not require repositioning and is a good choice for those who prefer nothing at their feet
- Non-Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Length	Tubing	Qty/Case
94010	WaterBoom® Floor Suction Strip	32"	Yes	10
94030	WaterBoom® Floor Suction Strip	32"	Yes	30

ViraSorb® Solidifier

- Evenly gels canister contents in seconds without stirring or agitating
- Non-toxic, Glutaraldehyde-free, and can be safely incinerated
- It does not revert back to a fluid state, helping to eliminate spillage and leakage, qualifying it as Solid Regulated Medical Waste
- Pour nozzles are included with multi-use bottles
- Non-Sterile and not made with natural rubber latex



Part No.	Description	Solidifies	Qty/Case
1000	ViraSorb® Solidifier – Single-Use Bottle	1,200cc	96
1001	ViraSorb® Solidifier – Single-Use Bottle	1,500cc	96
1002	ViraSorb® Solidifier – Single-Use Bottle	2,000cc	36
1003	ViraSorb® Solidifier – Single-Use Bottle	3,000cc	36
1005	ViraSorb® Solidifier – Single-Use Pouch	500cc	200
1010	ViraSorb® Solidifier – Multi-Use Bottle	10,000cc	12
1020	ViraSorb® Solidifier – Multi-Use Bottle	20,000cc	12

Richard-Allan®, Securline®, and WriteSite® Surgical Marking Pens

- Intended to provide an accurate and reliable way to mark the incision site prior to surgery
- Contain gentian violet ink
- A variety of barrel and tip options are available depending on preference
- Intended for single-use and not made with natural rubber latex

Richard-Allan® Classic Regular Tip Skin Marker

Part No.	Description	Sterility	Qty
2550	Marker, Ruler, and Label Set	Sterile	12/Box
2551	Marker	Sterile	12/Box
2650	Marker, Ruler, and Label Set	Sterile	50/Box
2651	Marker	Sterile	50/Box

Richard-Allan® Classic Fine Tip Skin Marker

Part No.	Description	Sterility	Qty
2530	Marker, Ruler, and Label Set	Sterile	12/Box
2531	Marker	Sterile	12/Box
2630	Marker, Ruler, and Label Set	Sterile	50/Box
2631	Marker	Sterile	50/Box

Richard-Allan® Modern Regular Tip Skin Marker

Part No.	Description	Sterility	Qty
2712	Marker, Ruler, and Label Set	Sterile	12/Box
2750	Marker, Ruler, and Label Set	Sterile	50/Box
2751	Marker	Sterile	50/Box

Richard-Allan® Taper Tip Skin Marker

Part No.	Description	Sterility	Qty
2730	Marker, Ruler, and Label Set	Sterile	50/Box
2731	Marker	Sterile	50/Box

Securline® Taper Tip Skin Marker

Part No.	Description	Sterility	Qty
1000-00-PDG	Marker and Ruler Set	Sterile	10/Box, 10 Boxes/Case
100L-00-PDG	Marker, Ruler, and Label Set	Sterile	100/Box
1002-00-PDG	Marker and Ruler Set	Sterile	100/Box
1003-00-PDG	Marker and Ruler Set	Non-Sterile	100/Box

WriteSite® Preoperative Surgical Site Marker

Part No.	Description	Sterility	Qty
2753BN	Marker Only	Non-Sterile	100/Box

WriteSite® Jr. Small Preoperative Skin Marker

Part No.	Description	Sterility	Qty
2900BN	Marker Only	Non-Sterile	100/Box

WriteSite® Dual Skin Marker


Part No.	Description	Sterility	Qty
2777	Marker Only	Sterile	50/Box
2778BN	Marker Only	Non-Sterile	100/Box

WriteSite® Multi-Ink Marker

Part No.	Description	Sterility	Qty
2754	Marker and 9 Label Set	Sterile	12/Box
2754BN	Marker and 6 Label Set	Non-Sterile	100/Box
2755	Marker Only	Sterile	12/Box
2755BN	Marker Only	Non-Sterile	100/Box

WriteSite®Plus Surgical Skin Markers

- Aspen Surgical is the only manufacturer to carry a line of surgical marking pens that include an ink specifically designed for use with ChloraPrep®, sold under the WriteSite® Plus brand
- Ink is also compatible with Duraprep® and Betadine® prep solutions
- Intended for single-use and not made with natural rubber latex



WriteSite® Plus Standard Size Surgical Skin Marker



Part No.	Description	Sterility	Qty
2700	Marker Only	Sterile	50/Box
2701	Marker, Ruler, and Label Set	Sterile	50/Box
2700BN	Marker Only	Non-Sterile	100/Box
2705BN	Marker Only	Non-Sterile	50/Jar

WriteSite® Plus Jr Size Surgical Skin Marker



Part No.	Description	Sterility	Qty
2715	Marker Only	Sterile	50/Box
2716	Marker, Ruler, and Label Set	Sterile	50/Box
2715BN	Marker Only	Non-Sterile	100/Jar

Securline® Laboratory Marker

Securline® Fine Tip Lab Marker

- Used to make permanent markings on plastic, metal, and porcelain labware
- Can be used on damp surfaces
- Available in four different colors
- Intended for single-use and not made with natural rubber latex



Part No.	Description	Color	Sterility	Qty/Box
1400-13-PDC	Marker Only	Blue	Non-Sterile	10
1400-15-PDC	Marker Only	Green	Non-Sterile	10
1400-16-PDC	Marker Only	Red	Non-Sterile	10
1400-20-PDC	Marker Only	Black	Non-Sterile	10

Securline® SuperFrost® Lab Marker

- Resistant to water and alcohol
- Recommended on frosted glass and plastic surfaces
- Available in two different colors
- Intended for single-use and not made with natural rubber latex



Part No.	Description	Color	Sterility	Qty/Box
1450-16-PDC	Marker Only	Red	Non-Sterile	10
1450-20-PDC	Marker Only	Black	Non-Sterile	10

Richard-Allan® Utility Markers

- Intended to provide a permanent method of identifying instruments, specimens, medication labels, and other objects on and off the sterile field
- Contains a black, waterproof ink that can withstand autoclave sterilization and can be used on most surfaces
- Available in several barrel and tip options and with or without labels
- Intended for single-use and not made with natural rubber latex

Richard-Allan® Regular Tip Utility Marker



Part No.	Description	Sterility	Qty/Box
2850	Marker and 9 Label Set	Sterile	50
2850BN	Marker and 6 Label Set	Non-Sterile	100
2851	Marker Only	Sterile	50
2851BN	Marker Only	Non-Sterile	100

Richard-Allan® Taper Tip Utility Marker



Part No.	Description	Sterility	Qty/Box
2830BN	Marker and 6 Label Set	Non-Sterile	100
2831	Marker Only	Sterile	50
2831BN	Marker Only	Non-Sterile	100

Select markers are available with Blank Labels, Rulers, or Medication Labels as noted.



9 Black and White Medication Labels



Blank Labels



Rulers

Medication Labels

- Intended to help easily identify and record medications in the O.R. for specific departments
- Color-coded to aid in easy identification, including medication name, strength, amount, and expiration date
- Labels are waterproof, smear-proof, and will not leave residue
- Twenty labels on each sheet
- Available sterile and non-sterile and with or without Richard-Allan® Utility Markers to help aid in convenient labeling
- Intended for single-use and not made with natural rubber latex



Medication Label Only

Part No.	Description	Sterility	Qty/Sheet
ML5840BN	Cath Lab	Non-Sterile	1000
ML5841BN	General Surgery 1	Non-Sterile	1000
ML5842BN	General Surgery 2	Non-Sterile	1000
ML5843BN	Cardiology	Non-Sterile	1000
ML5844BN	Cardiac	Non-Sterile	1000
ML5845BN	Neurology	Non-Sterile	1000
ML5847BN	Eyes	Non-Sterile	1000
ML5848BN	Ophthalmology	Non-Sterile	1000
ML5849BN	Pain	Non-Sterile	1000
ML5850BN	Interventional Radiology	Non-Sterile	1000
ML5851BN	Surgery Center	Non-Sterile	1000
ML5852BN	ENT	Non-Sterile	1000

Medication Label with Richard-Allan® Utility Marker

Part No.	Description	Sterility	Qty/Sheet
ML5832	Standard	Sterile	50
ML5840	Cath Lab	Sterile	50
ML5841	General Surgery 1	Sterile	50
ML5842	General Surgery 2	Sterile	50
ML5843	Cardiology	Sterile	50
ML5844	Cardiac	Sterile	50
ML5845	Neurology	Sterile	50
ML5847	Eyes	Sterile	50
ML5848	Ophthalmology	Sterile	50
ML5849	Pain	Sterile	50
ML5850	Interventional Radiology	Sterile	50
ML5851	Surgery Center	Sterile	50
ML5852	ENT	Sterile	50



Syringe and Scalpel Handling Systems

- Used to provide safe and easy handling and storage for syringes and scalpels
- Strong adhesive backed block that firmly grips to drapes
- Sterile, intended for single-use, and not made of natural rubber latex

Foam Syringe Stand



Part No.	Block Size	Width mm	Width Inches	Qty/Box
Sterile				
111001PBX	Single	31.8	1.25	5
111002PBX	Double	44.5	1.75	5
111003PBX	Triple	69.9	2.75	5

Foam Scalpel Holder



Part No.	Block Size	Width mm	Width Inches	Qty/Box
Sterile				
112001PBX	One Size	31.8	1.25	5

AspenBasin® Fluid Waste Container

- Used for general disposal of waste fluids in the Catheterization Lab
- Helps reduce exposure to contaminated fluid and bloodborne pathogens
- Grooved disposal area allows for convenient placement of syringe during disposal
- Includes high fluid capacity exceeding 500ml
- Secure lid design allows for easy disposal
- Sterile, intended for single-use, and not made of natural rubber latex



Part No.	Description	Qty
AB100	AspenBasin® Fluid Waste Container	10

Bard-Parker® Scalpel Holder

- Provides a safe placement of sharps between handling
- Helps prevent accidental sharps injuries, while keeping scalpels ready
- Includes slots for three scalpels to be held in a safe position
- Includes an adhesive back to ensure a safe, secure, and convenient placement in the sterile field
- Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Qty/Box
ST-100	Bard-Parker® Scalpel Holder	100

Bard-Parker® Hands-Free Transfer Tray

- Provides greater security than bowls offer on uneven surfaces
- Designed with many carrying points and shields your fingers providing ultimate protection
- Designed to hold most commonly used scalpel and blade sizes
- Strong, durable and light weight to provide easy handling
- Bright red for high visibility to sharp instruments
- Optimized for easy hand-free transfer of sharp instruments
- Helps minimize exposure to sharp instruments
- Sterile, intended for single-use, and not made with natural rubber latex




Part No.	Description	Qty/Case
ST-110	Bard-Parker® Hands-Free Transfer Tray	24

Bard-Parker™ Basin and Tray Liner Transfer Systems

- Provides a highly visible non-slip surface to help promote O.R. safety and instrument care
- Basin liners are available for large and medium sized kidney basins
- Tray liners are available for 10" x 6" x 2" and 8" x 5" x 2" trays
- Available in formable and non-formable versions
- Formable versions allow for personalized configuration
- Reusable
- Non-sterile and not made with natural rubber latex

NEW!



Part No.	Description	Size	Qty/Box
ST-101NS	Basin Liner Transfer System		5
ST-102NS	Tray Liner Transfer System		5
ST-104NS	Basin Liner Transfer System		5
ST-151NS	Basin Liner Transfer System		5
ST-154NS	Basin Liner Transfer System		5

Bard-Parker™ Surgical Pause Mat and Transfer Zone Mat

- Helps provide a highly visible non-slip surface to promote O.R. safety and instrument care
- Helps support the Joint Commission's Universal Protocol for Surgical Time Out
- Reusable
- Non-sterile and not made with natural rubber latex

Transfer Zone Mat

- Available in 8" x 11" and 8" x 14" sizes
- Offered in a formable version for personalized configuration and sharps handling

NEW!



Part No.	Description	Size	Color	Qty/Box
ST-202NS	Transfer Zone Mat	8" x 11"	Yellow	5
ST-203NS	Transfer Zone Mat	8" x 14"	Yellow	5
ST-302NS	Transfer Zone Mat	8" x 11"	Yellow	10
ST-303NS	Transfer Zone Mat	8" x 14"	Yellow	10

NEW!



Surgical Pause Mat

- Available in a 8" x 11" size
- Use Surgical Pause side as a reminder and use the opposite side as a neutral zone
- Available in 2 color options

Part No.	Description	Size	Color	Qty/Box
ST-304NS	Surgical Pause Mat	8" x 11"	Red/Green	10
ST-305NS	Surgical Pause Mat	8" x 11"	Red/Yellow	10

ErgoSupport™ Patient Transfer Board

The ErgoSupport™ Patient Transfer Board is a lightweight device designed to make lateral patient transfer easier for O.R. staff, and help make it more comfortable for the patient. It's designed to be more comfortable for patients than traditional metal roller boards because of its soft foam core. It is also lightweight, for easy transportability and storage. In addition, Aspen Surgical patient transfer boards are available in three sizes to accommodate different patient transfer needs throughout the hospital.

- Ideal for use with:
- Transferring patients from O.R. tables to stretchers
  - Transferring patients from stretchers to hospital beds
  - Comfortable – soft core with no hard metal rollers
  - Lightweight – no lifting required, even with your heavy patients<sup>1</sup>
  - Simple – uses same or similar transfer technique as those used with metal roller boards
  - Fowler Position Transfer – long transfer board allows for transfer in seated position
  - Foldable – long transfer board folds for easy storage and set-up
  - Strap – carry this lightweight device on your shoulder and keep your hands free
  - Resilient and sturdy – durable material protects against end cap breakage

Part No.	Description	Sterility	Qty/Box
93000	Long, Wide, 71.5" x 20.5"	Non-Sterile	1
93001	Short, Narrow, 40" x 16.5"	Non-Sterile	1
93002	Short, Wide, 32" x 20.5"	Non-Sterile	1



NEW!







Richard-Allan® Surgical Needles

Aspen Surgical carries many styles and sizes of Richard-Allan® Free-Eyed Disposable Needles used for suturing. This includes needles used in an operating room setting. Richard-Allan® Needles are known for their strength and are hand-crafted to help guard against imperfections. This makes Richard-Allan® Surgical Needles a great choice for specialty procedures.


Please note the needles drawn are not to scale.

Surgical Needle Reference Guide

2050 – Reference Number

Description – Kidney Blunt ½ Circle Taper

Eye/Point Diagram




Wire Diameter

205000 0.056 "

Needle Length

2.244 "



Needle

Straight Eye		Reverse Cutting	Cutting	Taper	Trocar	Blunt
Regular Surgeon's			2134 VL-M	2146 L	None	2051 H
				2147 VL		213413M
				2200 VL		

⅔ Circle Eye		Reverse Cutting	Cutting	Taper	Trocar	Blunt
Regular Surgeon's		2090 VL-H	2441 H	2168 H	2169 H	None
		2124 VL		2195 M		
		2145 M				
		2153 VL-L				
		2158 VL				

½ Circle Eye		Reverse Cutting	Cutting	Taper	Trocar	Blunt
Regular Surgeon's		2091 VL-M	None	2092 VL	2170 H	2050 H
		2123 VL		2095 L-H	2179 M-H	
		2127 L-H		2116 M	2426 M-H	
		2128 L-H		2125 VL	2701 L	
		2135 M-H		2140 L-H		
		2154 VL-L		2150 VL		
		2157 VL		2151 L		
		2164 M-VH		2159 L-M		
		2432 L-M		2167 H		
				2180 M		
				2181 L-M		
				2222 M		
				2442 VL-L		
				2500 VL		
				2706 L		
				2707 L		

Round	None			2126 VL-H	None	None
				2223	VH	


Spring	2156 VL	2000 VL	2129 VL	None	None
			2002 M		

Scale						
VL	Very Light Wire	L	Light Wire	M	Medium Wire	H
0.021"		0.027"		0.037"		0.050"
0.024"		0.030"		0.040"		0.056"
		0.034"		0.043"		

Richard-Allan® Surgical Needles


2000


French Spring Eye ½ Circle Cutting



200016 0.024"


1.024"





200018 0.024"


1.024"



Cutting Spring Eye


2002


French Spring Eye Kelly Intestinal ½ Circle Taper



200202 0.035"


1.417"





200204 0.031"


1.102"



Taper Double Spring Eye


2051

Kidney Blunt Straight



205100 0.056"

2.244"




Blunt Regular




Not to scale

Note: not to scale.






Reverse Regular




Reverse Regular



2092

Adson Dura ½ Circle Taper Point



209201 0.024"  
0.748"



209203 0.024"  
0.512"




209202 0.024"  
0.630"



2095

Ferguson (Regular Surgeon’s Eye) ½ Circle Taper Point



209501 0.056"  
3.110"



209510 0.040"  
1.693"



209503 0.056"  
2.677"



209511 0.037"  
1.535"



209504 0.050"  
2.480"



209512 0.037"  
1.535"



209505 0.050"  
2.323"



209513 0.037"  
1.378"



209506 0.050"  
2.323"



209515 0.034"  
1.220"



209507 0.043"  
2.087"



209516 0.03"  
1.220"



209508 0.043"  
2.087"



209517 0.030"  
1.063"



209509 0.040"  
1.890"




209518 0.03"  
1.063"



2116

Butler’s Cervix ½ Circle Taper Point



211601 0.037"  
0.984"




211602 0.037"  
1.378"



Taper Regular

2125

Eye Needle ½ Circle Taper Point



212501 0.024"  
1.063"



212503 0.021"  
0.748"



212502 0.024"  
0.866"



212504 0.021"  
0.630"




Taper Regular



2126

Ferguson’s (Round Eye) ½ Circle Taper Point



Taper Round


212601 0.056"  
3.110"




212602 0.056"  
2.953"




212603 0.050"  
2.638"




212604 0.050"  
2.480"




212605 0.043"  
2.362"




212606 0.043"  
2.205"



212607 0.040"  
1.890"



212608 0.037"  
1.732"



212609 0.037"  
1.496"



212610 0.034"  
1.299"



212611 0.034"  
1.299"



212612 0.030"  
1.024"



212613 0.027"  
0.866"



212614 0.024"  
0.827"



212615 0.024"  
0.669"



212616 0.024"  
0.669"




212622 0.024"  
1.102"



2127

Fistula Light ½ Circle Reverse Cutting



Reverse Regular

212703 0.043"  
2.087"



212704 0.037"  
1.693"



212705 0.034"  
1.496"



212706 0.034"  
1.220"




212708 0.027"  
0.827"



2128

Fistula Heavy ½ Circle Reverse Cutting




Reverse Regular


212801 0.056"  
2.362"



212802 0.056"  
2.126"



212803 0.050"  
1.850"



212804 0.043"  
1.654"



212805 0.037"  
1.378"




212806 0.034"  
1.102"



2129

French Spring Eye Intestinal ½ Circle Taper Point



Taper Double Spring Eye

212901 0.012"  
1.181"



212902 0.024"  
1.181"



212903 0.024"  
0.984"



212904 0.024"  
0.866"



212905 0.012"  
0.709"




212906 0.012"  
0.591"



2134

Keith’s Abdominal Straight Cutting



Cutting Regular

213401 0.040"  
4.016"



213402 0.040"  
3.465"



213404 0.034"  
2.992"



213405 0.034"  
2.756"



213406 0.034"  
2.480"



213407 0.030"  
2.244"



213408 0.030"  
2.008"



213409 0.030"  
1.732"



213410 0.027"  
1.496"




213412 0.024"  
1.024"





Note: 213413 = Blunt



2135

Kelly's Cervix ½ Circle Reverse Cutting




Reverse Regular

<b>213501</b>	<b>0.056"</b>	
	<b>2.126"</b>	
	<b>0.050"</b>	
	<b>1.693"</b>	



<b>213503</b>	<b>0.043"</b>	
	<b>1.654"</b>	
<b>213504</b>	<b>0.037"</b>	
	<b>1.378"</b>	


2140

Kelly's Intestinal ½ Circle Taper Point




Taper Regular

<b>214001</b>	<b>0.037"</b>	
	<b>1.535"</b>	
<b>214002</b>	<b>0.037"</b>	
	<b>1.378"</b>	




<b>214004</b>	<b>0.034"</b>	
	<b>1.102"</b>	

2145

King's ⅔ Circle Reverse Cutting




Reverse Regular






<b>214506</b>	<b>0.050"</b>	
	<b>4.016"</b>	
<b>214505</b>	<b>0.043"</b>	
	<b>3.740"</b>	
<b>214504</b>	<b>0.040"</b>	
	<b>3.465"</b>	

2146

Mayo Intestinal Straight Taper Point




Taper Regular






<b>214601</b>	<b>0.034"</b>	
	<b>2.992"</b>	
<b>214602</b>	<b>0.034"</b>	
	<b>2.756"</b>	
<b>214603</b>	<b>0.030"</b>	
	<b>2.008"</b>	
<b>214604</b>	<b>0.030"</b>	
	<b>1.732"</b>	
<b>214605</b>	<b>0.027"</b>	
	<b>1.220"</b>	

2147

Fine Intestinal Straight Taper Point




Taper Regular

<b>214701</b>	<b>0.024"</b>	
	<b>2.244"</b>	
<b>214702</b>	<b>0.024"</b>	
	<b>1.969"</b>	
<b>214703</b>	<b>0.024"</b>	
	<b>1.732"</b>	
<b>214704</b>	<b>0.024"</b>	
	<b>1.496"</b>	
<b>214706</b>	<b>0.024"</b>	
	<b>1.102"</b>	

2151

Mayo Intestinal ½ Circle Taper Point




Taper Regular

<b>215102</b>	<b>0.027"</b>		<b>215105</b>	<b>0.027"</b>	
	<b>1.732"</b>			<b>1.220"</b>	
			<b>215106</b>	<b>0.027"</b>	
				<b>1.063"</b>	



2153



Reverse Regular

Lane's Cleft Palate 3/8 Circle Reverse Cutting

215303 0.024"  
0.748"



215305 0.024"  
0.512"




215304 0.024"  
0.630"



215306 0.024"  
0.394"



2154



Reverse Regular

Lane's Cleft Palate 1/2 Circle Reverse Cutting

215402 0.027"  
0.827"



215405 0.024"  
0.512"



215403 0.024"  
0.748"




215406 0.024"  
0.394"



215404 0.024"  
0.630"



2156



Reverse Double Spring Eye

Lane's 1/2 Circle Reverse Cutting French Spring Eye

215601 0.012"  
0.984"



215603 0.012"  
0.709"




215602 0.012"  
0.827"



215604 0.012"  
0.591"



2157



Lane's 1/2 Circle Reverse Cutting

215701 0.024"  
1.024"



215703 0.024"  
0.748"




215702 0.024"  
0.866"



215704 0.024"  
0.630"



2158




Reverse Regular

Lane's 3/8 Circle Reverse Cutting

215803 0.024"  
0.748"



2159



Taper Regular

Ligating 1/2 Circle Taper Point


215902 0.034"  
1.496"



215903 0.034"  
1.339"



2164



Reverse Regular

Martin's Uterine 1/2 Circle Reverse Cutting

216401 0.064"  
2.953"



216404 0.050"  
1.693"



216402 0.064"  
2.480"



216405 0.050"  
1.496"



216403 0.056"  
1.969"




216406 0.043"  
1.102"



216407 0.037"  
0.984"



2167



Taper Regular

Mayo Catgut 1/2 Circle Taper Point

216700 0.056"  
2.244"



216735 0.050"  
1.417"



216701 0.056"  
1.969"



216704 0.050"  
1.260"



216702 0.050"  
1.732"



216705 0.050"  
1.102"



216703 0.050"  
1.535"




216707 0.050"  
0.906"



2168

Mayo Catgut 3/8 Circle Taper Point (Heavy)



216801 0.056"  
1.969"


216803 0.050"  
1.535"

216802 0.050"  
1.732"

216804 0.050"  
1.260"

2169

Mayo Catgut 3/8 Circle Trocar Point



216901 0.056"  
1.969"

216904 0.050"  
1.260"

216902 0.050"  
1.732"


216905 0.050"  
1.102"

216903 0.050"  
1.535"

216906 0.050"  
0.984"

2170

Mayo Catgut 1/2 Circle Trocar Point



217001 0.056"  
1.969"

217004 0.05"  
1.260"

217002 0.05"  
1.732"


217005 0.05"  
1.102"

217003 0.05"  
1.535"

217006 0.05"  
0.984"

2179

Murphy's Trachelorrhaphy 1/2 Circle (Light Trocar)



217901 0.050"  
1.260"


217903 0.043"  
0.984"

217902 0.050"  
1.102"

217904 0.037"  
0.748"

2180


Murphy's Hemostatic 1/2 Circle Taper Point



218000 0.037"  
0.984"

2181

Murphy's Intestinal 1/2 Circle Taper Point



218101 0.037"  
1.693"


218103 0.034"  
1.339"

218102 0.037"  
1.535"

218104 0.034"  
1.220"

2195


Davis Tonsil



219500 0.043"  
0.984"

2200

Bunnels Straight



220004 0.024"  
2.480"


220003 0.021"  
2.244"

220002 0.021"  
1.969"



2222

½ Circle Improved Mayo Catgut Taper Point (Light)



222201 0.043"  
1.969"



222203 0.043"  
1.535"



222204 0.043"  
1.260"



222202 0.043"  
1.693"




222205 0.043"  
1.102"



2223

Gallies Fascia ½ Circle Taper



Taper Round

222301 0.080"  
2.520"




222302 0.080"  
2.008"



2426

Mayo Catgut ½ Circle Trocar



Trocar Regular

242602 0.050"  
1.969"



242605 0.043"  
1.378"



242604 0.050"  
1.535"




242606 0.040"  
1.102"



2432

Fistula ½ Circle Reverse Cutting



Reverse Regular

243204 0.043"  
1.654"



243207 0.037"  
0.984"



243205 0.043"  
1.378"



243208 0.034"  
0.866"



243206 0.040"  
1.260"




243209 0.030"  
0.748"



2441

King Cutting ¾ Circle




Cutting Regular

244101 0.050"  
3.465"



2442

Ferguson's ½ Circle Taper Point



Taper Regular

244216 0.027"  
1.063"



244220 0.024"  
0.630"



244218 0.024"  
0.748"




244222 0.024"  
0.512"



2500

½ Circle Taper




Taper Regular

250000 0.024"  
1.063"



2706

Ferguson ½ Circle Taper




Taper Regular

270616 0.030"  
1.063"



2707

Ferguson ½ Circle Taper



Taper Regular

270718 0.027"  
0.787"



- Saber® 100ml Infection Control  
PVC Bulb Evacuator

  - Advanced vacuum level to help prevent occlusion and premature clotting off
  - Large lumen “maxi-flow” anti-reflux valve on the in-flow port
  - Provides exceptional in-flow, even when clots are present
  - Large lumen “maxi-flow” anti-reflux valve on the bottom emptying port
  - Emptying port can be connected to auxiliary suction
  - Easy to empty with reduced chance of bacterial migration into the bulb evacuator
  - Sterile, intended for single-use, and not made with natural rubber latex



- 3C™ Collection System

  - Maintains a completely closed system to help reduce chance of infection
  - Keeps the closed system intact, even when emptying bulb and replacing the collection system
  - Provides added fluid collection capacity with 150ml or 500ml 3C Collection Systems
  - Anti-reflux valve helps reduce the chance of exposure to blood borne pathogens during disposal
  - Sterile, intended for single-use, and not made with natural rubber latex



Bulb Evacuator and Collection Systems		
Part No.	Description	Qty/Case
310100	Saber® 100ml Infection Control PVC Bulb Evacuator with 150ml 3C™ Collection System	10
310101	Saber® 100ml Infection Control PVC Bulb Evacuator	10

- S-Vac® 100ml PVC Bulb Evacuator

  - Advanced vacuum level to help prevent occlusion and premature “clotting off”
  - Large lumen “maxi-flow” anti-reflux valve on the in-flow port
  - Provides exceptional in-flow, even when clots are present
  - Unique cap/plug emptying port plug closure with external skirt
  - Helps prevent touch contamination of the port
  - Emptying port can be connected to auxiliary suction
  - Appropriate option when higher level suction is needed
  - Sterile, intended for single-use, and not made with natural rubber latex



Bulb Evacuator		
Part No.	Description	Qty/Case
312101	S-VAC® 100ml PVC Bulb Evacuator	10

- S-Vac® Silicone Bulb Evacuator


  - Provides lower vacuum level than the Saber® 100 or the S-Vac® PVC bulbs and is preferred by some physician specialties
  - Sterile, intended for single-use, and not made with natural rubber latex




Bulb Evacuator		
Part No.	Description	Qty/Case
322160	S-VAC® Silicone 100ml Bulb Evacuator	10

- S-Vac® Fluted Channel Drains


  - Designed to provide advanced drainage and increased patient comfort
  - Fluted channel design helps minimize the chance of tissue in-growth creating the blockages that can cause premature drain removal when using standard drains
  - 10mm Flat Fluted Drain provides 6-channels for optimal drainage, compared to the industry standard 4-channel design
  - Soft, flat drain hubs offers enhanced patient comfort
  - Exceptional flow rates provides improved drainage, helping reduce the chance of seroma formation
  - Available with or without sharp stainless steel trocars
  - Radiopaque
  - Sterile, intended for single-use, and not made with natural rubber latex



6 Channel Fluted Flat Drain



4 Channel Fluted Flat Drain



4 Channel Fluted Round Drain

Fluted Flat Drain		
Part No.	Description	Qty/Case
322211	7mm Flat, Full Fluted, 4 Channel Drain	10
322213	10mm Flat, ¾ Fluted, 6 Channel Drain	10
322214	10mm Flat, Full Fluted, 6 Channel Drain	10

Fluted Flat Drain with Trocar		
Part No.	Description	Qty/Case
322212	7mm Flat, Full Fluted, 4 Channel Drain, Trocar	10
322215	10mm Flat, Full Fluted, 6 Channel Drain, Trocar	10
322216	7mm Flat, ¾ Fluted, 4 Channel Drain, Trocar	10
322217	10mm Flat, ¾ Fluted, 6 Channel Drain, Trocar	10

Fluted Round Drain		
Part No.	Description	Qty/Case
332186	10 FR Round, Full Fluted, 4 Channel Drain	10
332188	15 FR Round, Full Fluted, 4 Channel Drain	10
332190	19 FR Round, Full Fluted, 4 Channel Drain	10

Fluted Round Drain with Trocar		
Part No.	Description	Qty/Case
332187	10 FR Round, Full Fluted, 4 Channel Drain, Trocar	10
332189	15 FR Round, Full Fluted, 4 Channel Drain, Trocar	10
332191	19 FR Round, Full Fluted, 4 Channel Drain, Trocar	10
332221	10 FR Round, ¾ Fluted, 4 Channel Drain, Trocar	10
332223	15 FR Round, ¾ Fluted, 4 Channel Drain, Trocar	10
332225	19 FR Round, ¾ Fluted, 4 Channel Drain, Trocar	10



6 Channel Flat



4 Channel Flat



4 Channel Round


Xtra Fluted Round Drain		
Part No.	Description	Qty/Case
322226	10 FR Round Xtra Fluted Drain	10
322227	10 FR Round Xtra Fluted Drain, With Trocar	10
322228	15 FR Round Xtra Fluted Drain	10
322229	15 FR Round Xtra Fluted Drain, With Trocar	10
322230	19 FR Round Xtra Fluted Drain	10
322231	19 FR Round Xtra Fluted Drain, With Trocar	10

Drain System-Fluted with 100ml Bulb Evacuator		
Part No.	Description	Qty/Case
332259	10 FR Drain System Round ¾ Fluted, with Trocar and 100ml S-Vac® Bulb Evacuator	10
332263	10 FR Drain System Round Full Fluted, with Trocar and 100ml S-Vac® Bulb Evacuator	10
332264	15 FR Drain System Round Full Fluted, with Trocar and 100ml S-Vac® Bulb Evacuator	10
332265	19 FR Drain System Round Full Fluted, with Trocar and 100ml S-Vac® Bulb Evacuator	10



Perforated Drains

- Discrete drain with our unique channel designed to help prevent collapsing or kinking
- Discrete drains contain internal ribs to help prevent collapse and occlusion of the drain
- 10mm discrete drain offers minimal transition from the drain hub to the drain tubing for ease of removal and added patient comfort
- Hubless styles also available
- Soft, silicone material
- Available with or without extra sharp stainless steel trocars
- Radiopaque
- Sterile, intended for single-use, and not made with natural rubber latex



Flat Perforated Drain

Round Perforated Drain

Discrete Flat Perforated Drain

Part No.	Description	Qty/Case
370004	Discrete Drain, Flat 7mm x 20cm, ¾ Perf	10
370005	Discrete Drain, Flat 10mm x 20cm, ¾ Perf	10
370006	Discrete Drain, Flat 7mm x 20cm, Full Perf	10
370007	Discrete Drain, Flat 10mm x 20cm, Full Perf	10

Round Perforated Drain

Part No.	Description	Qty/Case
370020	Discrete Drain, Round 7 FR, End-Perf	10
370021	Discrete Drain, Round 10 FR, End-Perf	10
370022	Discrete Drain, Round 15 FR, End-Perf	10
370023	Discrete Drain, Round 19 FR, End-Perf	10

Discrete Flat Perforated Drain With Trocar

Part No.	Description	Qty/Case
370000	Discrete Drain, Flat 7mm x 20cm, ¾ Perf, Trocar	10
370001	Discrete Drain, Flat 10mm x 20cm, ¾ Perf, Trocar	10
370002	Discrete Drain, Flat 7mm x 20cm, Full Perf, Trocar	10
370003	Discrete Drain, Flat 10mm x 20cm, Full Perf, Trocar	10

Round Perforated Drain With Trocar

Part No.	Description	Qty/Case
370016	Discrete Drain, Round 7 FR, End-Perf, Trocar	10
370017	Discrete Drain, Round 10 FR, End-Perf, Trocar	10
370018	Discrete Drain, Round 15 FR, End-Perf, Trocar	10
370019	Discrete Drain, Round 19 FR, End-Perf, Trocar	10

Hubless Flat Perforated Drain

Part No.	Description	Qty/Case
360004	Hubless Drain, Flat 7mm x 20cm, ¾ Perf	10
360005	Hubless Drain, Flat 10mm x 20cm, ¾ Perf	10
360006	Hubless Drain, Flat 7mm x 20cm, Full Perf	10
360007	Hubless Drain, Flat 10mm x 20cm, Full Perf	10

Drain Systems with S-Vac® 100ml Bulb Evacuator

Part No.	Description	Qty/Case
340001	10mm Drain System, Flat ¾ Perforated, with Trocar and 100ml S-Vac® Bulb Evacuator	10
340002	7mm Drain System, Flat Full Perforated, with Trocar and 100ml S-Vac® Bulb Evacuator	10
340003	10mm Drain System, Flat Full Perforated, with Trocar and 100ml S-Vac® Bulb Evacuator	10
340004	7mm Drain System, Flat ¾ Perforated, with 100ml S-Vac® Bulb Evacuator	10
340005	10mm Drain System, Flat ¾ Perforated, with 100ml S-Vac® Bulb Evacuator	10
340006	7mm Drain System, Flat Full Perforated, with 100ml S-Vac® Bulb Evacuator	10
340007	10mm Drain System, Flat Full Perforated, with 100ml S-Vac® Bulb Evacuator	10
340018	15 FR Drain System, Round End-Perforated, with Trocar and 100 ml S-Vac® Bulb Evacuator	10

Flat Perforated Drain With Trocar

Part No.	Description	Qty/Case
360000	Hubless Drain, Flat 7mm x 20cm, ¾ Perf, Trocar	10
360002	Hubless Drain, Flat 7mm x 20cm, Full Perf, Trocar	10
360003	Hubless Drain, Flat 10mm x 20cm, Full Perf, Trocar	10

Secure-It™ Nasal Tube Holder

- Helps reduce the time spent securing nasogastric tubing and oxygen cannulas
- Increases patient comfort by holding tubing and cannulas at an optimum angle every time, minimizing the chance of tissue necrosis
- Flexible and breathable material will stay in place, even in the presence of moisture
- Non-Sterile, intended for single-use, and not made with natural rubber latex

Secure-It™ Nasal Tube Holder

Part No.	Description	Qty/Box
630110xbx	Large Adult	100

Transparent Dressings

- Intended for wound care applications, securing objects to the skin such as catheters, and as a secondary dressing
- Used to help protect the wound site from possible bacterial and viral infection
- Breathable to ensure the wound is kept moist
- Patented thumb tab design to make application easier
- Sterile, intended for single-use, and not made with natural rubber latex



Transparent Dressings

Part No.	Description	Qty/Box
TDS-24	2.375" x 2.75"	400 (100/box, 4 boxes/case)
TDS-26	4" x 4.5"	200 (50/box, 4 boxes/case)

## Sterion® Filament Reinforced Skin Closures


- Reinforcing filaments offer added strength for large skin surface closures
- Filament Reinforced Skin Closure style is accepted as the standard by most healthcare professionals
- Breathable material designed to allow moisture transmission, to help reduce the chance of skin maceration
- Less painful to apply than suturing or stapling providing more comfort to the patient
- No adhesive enhancers required
- Strips are curl-resistant
- Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	mm	inches	Strips/Pouch	Pouches/Box
751040pbx	Filament Reinforced Skin Closures	3 x 75	⅜ x 3	5	50
751041pbx	Filament Reinforced Skin Closures	6 x 75	¼ x 3	3	50
751042pbx	Filament Reinforced Skin Closures	6 x 38	¼ x 1½	6	50
751046pbx	Filament Reinforced Skin Closures	6 x 100	¼ x 4	10	50
751047pbx	Filament Reinforced Skin Closures	12 x 100	½ x 4	6	50
751048pbx	Filament Reinforced Skin Closures	25.4 x 127	1 x 5	4	25

## Sterion® Flexible Skin Closures

- Non-woven flexible backing helps reduce the chance of skin shear that can occur when edema or flexing puts excessive pressure on the skin closure ends
- Less painful to apply than suturing or stapling, providing more comfort to the patient
- No adhesive enhancers required
- Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	mm	inches	Strips/Pouch	Pouches/Box
753040pbx	Blend Tone Skin Closures	3 x 75	⅜ x 3	5	50
753041pbx	Blend Tone Skin Closures	6 x 75	¼ x 3	3	50
753042pbx	Blend Tone Skin Closures	6 x 38	¼ x 1½	6	50
753046pbx	Blend Tone Skin Closures	6 x 100	¼ x 4	10	50
753047pbx	Blend Tone Skin Closures	12 x 100	½ x 4	6	50
753048pbx	Blend Tone Skin Closures	25 x 125	1 x 5	4	25

### Notes:

[illegible]



Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world, Hill-Rom's people, products, and programs work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2016 Aspen Surgical, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-1010427 rev 5 21-MAR-2016 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

**USA 616-698-7100**

**Canada 888-364-7004**

**[www.aspensurgical.com](http://www.aspensurgical.com)**