

# SINGLE-USE CATALOG



Per Procedure Kits



Instrumentation



Corneal Transplant



Cataract Procedures



Vitreoretinal Procedures



Intraocular Lenses



Viscoelastic



Accessories



Suite 12, Bunkilla Plaza,  
Bracetown Business Park,  
Clonee, Dublin 15,  
D15XC64, Republic of Ireland.

Tel: + 353 1 806 8600  
Direct: + 353 1 834 4340  
Email: [info@tpwhelehan.ie](mailto:info@tpwhelehan.ie)  
Web: [www.whelehansurgical.ie](http://www.whelehansurgical.ie)

For Ordering, please contact:

Uniphar Hospital Services,  
4045 Kingswood Road,  
Citywest Business Park, Dublin 24.  
E-Mail: [Unipharhospitalservices@uniphar.ie](mailto:Unipharhospitalservices@uniphar.ie)

Tel: + 353 1 468 8456  
Fax: + 353 1 466 0782

# Table of Contents



## Per Procedure Kits

Select Pack	3
Custom Packs	4



## Instrumentation

Knives	5
Cannulas	8
Cataract & Vitreoretinal Instruments	10
I/A Handpieces	11



## Corneal Transplant

Corneal Storage	12
Trephines	13



## Cataract Procedures

1.8 mm MICS™ Procedure	14
2.2 mm Small Incision Procedure	15
2.8 mm Standard Procedure	16
Other Cataract Products	17



## Vitreoretinal Procedures

23 Gauge Procedure	18
25 Gauge Procedure	20
20 Gauge Procedure	22
Other Retina	24
Bipolar	25



## Intraocular Lenses

Crystalens Accommodating IOL	26
Akreos Monofocal IOL	27
SofPort IOL	28
PMMA Lenses	29



## Viscoelastics

Viscoelastics	30
---------------	----



## Accessories

Eye Shields	31
Corneal Shields	32
Miscellaneous	32





## Simplify Your Procedure Packs with the Select Pack System

Select Pack is a complete kit of components designed to streamline and optimize your cataract or vitreoretinal procedure. Combine this kit of the most popular procedural components with your choice of blades and cannulas to create a complete procedure kit.



### Select Pak Cataract



STERILE



Each Select Pack will come with standard components plus your choice of drape. Choose the drape you want in your Select Pack from the list below:

	Select Pack Catalog #	
	Incise	Oval Aperature
<b>FULL BODY</b>		
Ophthalmic Surgical Drape	SPK01	SPK01A
2 fluid pouches		
Breather Bar		
SMS Non-Woven Material		
71 x 100 inches		
<b>3/4 LENGTH</b>		
Ophthalmic Surgical Drape	SPK02	SPK02A
2 fluid pouches		
Breather Bar		
SMS Non-Woven Material		
55 x 71 inches		
<b>MID-SIZE</b>		
Ophthalmic Surgical Drape	SPK03	SPK03A
1 fluid pouch		
SMS Non-Woven Material		
48 x 55 inches		
Ophthalmic Surgical Drape	SPK04	SPK04A
1 fluid pouch		
Embossed Plastic Material		
48 x 51 inches		



### Select Pak Vit



STERILE



#### VITREORETINAL

Vitreoretinal Pack

Catalog #	Description
9/box	
MEY23VTBLA	30 Cotton Tip Wood Applicators, 3"
	3 Steri-Strip™ Reinforced Wound Closures
	10 Eye Spears
	3 Tube Holders with Velcro
	1 Instrument Pouch, 8 Pocket
	1 Rubberband
	3 Gauze Sponges, 4" X 4"
	1 Tie Twist
	1 Instrument Wipe
	6 Cotton Towels
	2 White Paper Towels
	1 Mayo Stand Cover
	1 Drape, Incise, 51" X 51"
	1 Drape, EENT Split, 76" X 124"
	2 Syringe 3ml, & 25g Needle
	1 Syringe, 1 ml
	1 Outer Wrap/Table Cover



# Custom Packs




## Build Your Custom Pack with these Simple Steps

STERILE



1. Determine your component list, including Storz® Ophthalmics Single-Use Instrumentation, B+L Ocudrapes and Laseredge disposable knives.
2. Contact B+L Customer Service at 800-338-2020 to initiate a quote.
3. A B+L Custom Packs Coordinator will review the request and work with you to ensure you are getting the best products at the best possible procedural price.
4. Once your Custom Designed Pack is finalized, begin receiving your inventory and spend less on handling and inventory control issues encountered when dealing with multiple suppliers.



## Sterile Drapes



STERILE



### FULL BODY

Ophthalmic Surgical Drape  
2 fluid pouches  
Breather Bar  
SMS Non-Woven Material  
71 x 100 inches

Incise  
DRA850A

Oval Aperature  
DRA851

### 3/4 LENGTH

Ophthalmic Surgical Drape  
2 fluid pouches  
Breather Bar  
SMS Non-Woven Material  
55 x 71 inches

DRA831

DRA830A

### MID-SIZE

Ophthalmic Surgical Drape  
1 fluid pouch  
SMS Non-Woven Material  
48 x 55 inches

DRA801

DRA800A

Ophthalmic Surgical Drape  
1 fluid pouch  
Embossed Plastic Material  
48 x 51 inches

DRA1061

DRA1059A

Ophthalmic Surgical Drape  
Embossed Plastic Material  
48 x 51 inches

DRA1060

### EQUIPMENT COVERS

Mayo Stand Cover  
Clear Plastic  
24 x 40 inches

50/box  
DRA7352



STERILE






## LASEREDGE® Knives



### E7590 Reusable Blade Tip Handle



## Guarded Knives for Safety




## Arrow Knives



STERILE



- 1.1 mm angled
- 0.8 mm angled



## Clear Corneal Knives



STERILE




### LASEREDGE®

- 1.3mm Angled Double Bevel
- 1.6mm Angled Double Bevel
- 2.75mm Angled Double Bevel
- 2.85mm Angled Double Bevel
- 2.9mm Angled Double Bevel
- 3.0mm Angled Double Bevel
- 3.2mm Angled Double Bevel

### GUARDED

- 2.75mm angled double bevel
- 2.85mm angled double bevel



## Crescent Knives



STERILE




### LASEREDGE®

- Straight
- Angled Bevel Up
- Angled Bevel Up Tip Only
- Angled Bevel Down

### GUARDED

- Angled Bevel Up



## Implant Knives



STERILE



- 3.2mm Angled Up
- 3.5mm Angled Up
- 4.0mm Angled Bevel Up
- 5.2mm Angled Bevel Up

Catalog #	Dimensions (mm)		
6/box	A	B	C
E7592	17.6	1.09	2.03
E7593	17.6	0.80	1.65


Catalog #	Dimensions (mm)		
6/box	A	B	C
E7594	12.6	1.3	1.45
E7595	12.6	1.6	1.7
E7527	14.4	2.75	3.1
E7528	14.4	2.85	3.1
E7587A	14.4	2.9	3.2
E7529	14.4	3.0	3.2
E7532	14.4	3.2	3.6
E7527G	18.5	2.75	3.10
E7528G	18.5	2.85	3.10

Catalog #	Dimensions (mm)		
6/box	A	B	C
E7500	12.7	2.1	3.7
E7510	12.7	2.1	3.7
E7510AT	12.7	2.1	3.7
E7510D	12.7	2.1	3.7
E7510G	16.8	2.1	3.7

Catalog #	Dimensions (mm)			
6/box	A	B	C	D
E7563	14.4	3.2	3.8	1.6
E7562	14.7	3.5	3.8	1.7
E7560A	15.2	4.0	3.4	1.2
E7561A	16.5	5.2	4.1	1.7



# Knives




## MVR Knives



STERILE



- 19 Gauge
- 20 Gauge
- 23 Gauge
- 25 Gauge




## Parallel Side Knives



STERILE



- 2.75mm Angled Double Bevel



## Slit Knives



STERILE



### LASEREDGE®

- 1.4mm Angled Bevel Up
- 1.5mm Angled Bevel Up
- 2.5mm Angled Bevel Up
- 2.65mm Angled Bevel Up
- 2.65mm Angled Bevel Up Tip Only
- 2.65mm Angled Double Bevel
- 2.75mm Angled Bevel Up
- 2.8mm Angled Double Bevel
- 2.85mm Angled Bevel Up
- 3.0mm Angled Bevel Up
- 3.0mm Angled Double Bevel
- 3.2mm Angled Bevel Up
- 3.2mm Straight

### GUARDED

- 2.5mm Angled Bevel Up
- 2.65mm Angled Bevel Up
- 2.65mm Angled Double Bevel
- 2.8mm Angled Double Bevel
- 2.85mm Angled Bevel Up
- 3.0mm Angled Bevel Up
- 3.2mm Angled Bevel Up


Catalog #	Dimensions (mm)		
6/box	A	B	C
E7519	31.8	1.5	3.1
E7520	31.8	1.2	2.4
E7523	31.8	0.65	1.1
E7525	31.8	0.5	1.0

Catalog #	Dimensions (mm)			
6/box	A	B	C	D
E7575	14.4	2.75	3.1	4.8

Catalog #	Dimensions (mm)		
6/box	A	B	C
E7598	13.3	1.4	1.4
E7597	13.3	1.5	1.5
E7550A	14.4	2.5	3.8
E7551A	14.6	2.6	3.9
E7551AT	14.6	2.6	3.9
E7551ADB	14.4	2.6	3.3
E7559A	14.4	2.75	3.2
E7552ADB	14.4	2.8	3.4
E7553A	14.1	2.8	3.3
E7556A	14.4	3.0	3.9
E7556ADB	14.4	3.0	3.7
E7557A	14.4	3.2	3.9
E7557	14.4	3.2	3.9

E7550AG	18.5	2.5	3.8
E7551AG	18.5	2.6	3.9
E7551ADBG	18.5	2.6	3.3
E7552ADBG	18.5	2.8	3.4
E7553AG	18.5	2.8	3.3
E7556AG	18.5	3.0	3.9
E7557AG	18.5	3.2	3.9





## Stab Knives



STERILE




### LASEREDGE®

- 15°
- 15° Tip Only
- 22.5°
- 30°
- 30° Tip Only
- 45°

### GUARDED

- 15°



## Trapezoidal Knives



STERILE



### LASEREDGE®

- 1.5-1.7mm Angled
- 1.5-2.0mm Angled
- 1.6-1.8mm Angled
- 2.0-2.2mm Angled
- 2.75-3.2mm Angled Bevel Up
- 2.85-3.2mm Angled Bevel Up

### GUARDED

- 1.6-1.8mm Angled
- 2.0-2.2mm Angled
- 2.85-3.2mm Angled Bevel Up

Catalog #	Dimensions (mm)		
6/box	A	B	C
E7515	11.2	1.6	5.9
E7515T	11.2	1.6	5.9
E7522	10.7	1.7	4.0
E7530	10.4	1.7	2.9
E7530T	10.4	1.7	2.9
E7545	8.9	1.9	1.9
E7515G	14.0	1.6	5.9

Catalog #	Dimensions (mm)				
6/box	A	B	C	D	E
E7599	14.4	N/A	2.0	2.0	N/A
E7596	14.4	2.0	1.8	1.5	
E7600	14.4	1.8	1.8	1.6	2.95
E7601	14.4	2.2	2.3	2.0	3.1
E7572	16.8	3.2	3.2	2.75	4.6
E7582	16.8	3.2	3.4	2.85	4.3
E7600G	18.5	1.8	1.8	1.6	2.95
E7601G	18.5	2.2	2.3	2.0	3.1
E7582G	18.5	3.2	3.4	2.85	4.3



# Cannulas



## Anesthesia Needles



STERILE



### RETROBULBAR

23 Gauge, Atkinson Style Bevel

25 Gauge, Atkinson Style Bevel

### PERIBULBAR

23 Gauge, Atkinson Style Bevel

25 Gauge, Atkinson Style Bevel

Catalog #	Dimensions (mm)	
5/box	Length (excluding hub)	
E7131		38.0
E7133		38.0
E7132		32.0
E7134		32.0



## Anterior Chamber Needles



STERILE



20 Gauge, Angled Tip

23 Gauge, Angled Tip

25 Gauge, Angled Tip

27 Gauge, Angled Tip

27 Gauge, Angled Tip

30 Gauge, Angled Tip

Catalog #	Dimensions (mm)	
5/box	Tip	Length (excluding hub)
E7101	7.0	24.5
E7102	7.0	24.5
E7103	9.0	20.5
E7104	5.0	20.5
E7106	5.0	21.0
E7105	5.0	21.0



## Cystotomes - Irrigating



STERILE



### CAPSULORHEXIS

25 Gauge, Formed Blunt Tip

27 Gauge, Formed Angled Tip

### FORMED

23 Gauge, Formed Curved Tip

25 Gauge, Formed Curved Tip

25 Gauge, Short Curved Tip

27 Gauge, Short Curved Tip

### OTHER

25 Gauge, Straight

27 Gauge, Curved Twisted Tip

Catalog #	Dimensions (mm)	
5/box	Tip	Length (excluding hub)
E7118		22.0
E7117		16.5
E7112	12.0	16.5
E7113	12.0	16.0
E7114	8.0	16.5
E7115	8.0	16.0
E7123		16.0
E7125 C		16.0



## Hydrodissection



STERILE



25 Gauge, Straight U-Tip Allows Irrigation/Aspiration at the 12 O'Clock Position

25 Gauge, Angled Flattened Tip for Insertion Under Capsular Flap

Catalog #	Dimensions (mm)	
5/box	Tip	Length (excluding hub)
E7154 S		19.0
E7155	7.0	21.0



## Lacrimal



STERILE



19 Gauge Shaft Tapered to 23 Gauge Tip, Straight

26 Gauge, Curved

Catalog #	Dimensions (mm)	
5/box	Length (excluding hub)	
E7184		37.5
E7187		31.0

# Cannulas



## Viscoelastic



STERILE



23 Gauge, Hub Facilitates Flow of  
Viscoelastic Material  
27 Gauge, Hub Facilitates Flow of  
Viscoelastic Material

Catalog #	Dimensions (mm)	
	Tip	Length (excluding hub)
5/box E7180	7.0	24.5
E7182	7.0	24.5



STERILE



## Vitreoretinal

### TRANSLOCATION

39 Gauge I.D., Flexible Tip

41 Gauge I.D., Flexible Tip

Catalog #	Description
5/box E7365	20 gauge shaft tapers to 39 gauge I.D. flexible tip. Used for hydrodissecting the retina for various macular hole procedures.
E7370	20 gauge shaft tapers to 41 gauge I.D. flexible tip. Used for hydrodissecting the retina for various macular hole procedures.



STERILE



### OTHER

23 Gauge, Soft Tip

25 Gauge, Soft Tip

10/box BL7823	Aspirating cannula with soft silicone tip.
BL7725	Aspirating cannula with soft silicone tip.



STERILE



23 Gauge Single Trocar/Cannula

20 Gauge Infusion

6/box CX9123	
MVS1038A	4.0 mm beveled tip and Luer Lock.



STERILE



## Miscellaneous

27 Gauge Silicone Capsule Polisher.  
Silicone tip provides squeegee effect.

25 Gauge Stromal Puncture Can-  
nula. For treatment of recurrent  
corneal erosion.

Catalog #	Dimensions (mm)	
	Tip	Length (excluding hub)
5/box E7153		22.0
E7185		22.0





# Cataract & Vitreoretinal Instruments



## Cataract Forceps



STERILE



Capsulorhexis Forceps

MICS™ Capsulorhexis Forceps

Bonn Forceps

Tying Forceps

McPherson Tying Forceps

Pierse Tying Forceps

Catalog #	Description
8/box SUE01	Straight jaws. Fine sharp tips.
SUE06	Fine angled tips for grasping delicate capsule during MICS 1.8 procedures.
SUE02	1x2 teeth 0.12 mm, 3.5 mm tying platform. For tissue manipulation and bulbus fixation
SUE03	3.5 mm tying platform. May be used for IOL manipulation and cartridge loading.
SUE04	Angled jaws, 4.0 mm tying platform. May be used for IOL cartridge loading.
SUE05	0.25 mm notched straight jaws. 4.0 mm tying platform. For tissue manipulation and bulbus fixation.



## Cataract Manipulators



STERILE



Koch Stop & Chop Manipulator

Sinskey Hook

Phaco Spatula

Catalog #	Description
8/box SUF01	Angled shaft with 1.5 mm tip. Designed to efficiently perform the "Stop and Chop" technique
SUF02	Angled shaft 10 mm. Tip length 0.5 mm. Diameter 0.18 mm. For lens dialing and manipulation.
SUF03	11.0 mm spatula length. For nucleus manipulation.



## Vitreoretinal Forceps



STERILE



### ENDGRIPPING

23 gauge

25 gauge

### ASYMMETRICAL PEELING

23 gauge

25 gauge

Catalog #	Description
6/box 65023S	2/box 65024S
65025S	65026S
SUI03	SUI04
SUI05	SUI06
Ultra-thin endgripping tip for enhanced visibility and grasp. Ergonomic handle consistent with our popular reusable design. Balanced design for optimal feel and control.	
Used for ILM and Epiretinal Membrane Peel. Use on macular puckers and cellophane maculopathy. Asymmetric jaw allows for excellent visualization of retinal tissue.	

# I/A Handpieces



## Silicone I/A Handpieces

12

STERILE



Silicone I/A Handpiece for the Stellaris® System

Silicone I/A Handpiece for the Millennium® and most competitive systems

Catalog #	Description
12/box 85910ST	Features a flexible 45° tip design for effective cortical removal. For use in a 2.2 - 2.8 mm incision. Smooth irrigation/aspiration ports eliminate sharp edges for reduced risk of capsule rupture. Silicone tip design facilitates cortex removal, capsule polishing, viscoelastic removal and IOL manipulation in the capsular bag. Semi-transparent silicone sleeve provides superior visualization.
85910S	



## I/A Handpieces

12

STERILE



**STELLARIS® SYSTEM**  
Curved I/A Handpiece

45° I/A Handpiece

90° I/A Handpiece

Curved Irrigation Handpiece

MICS™ 1.8 mm Curved I/A Handpiece with Sleeve

MICS™ 1.8 mm 45° I/A Handpiece with Sleeve

MICS™ 1.8 mm Curved Irrigation Handpiece

Nagahara Irrigating Chopper

### MILLENNIUM® & OTHER SYSTEMS

Curved Aspiration Handpiece

MICS™ 1.8 mm Curved Aspiration Handpiece

Catalog #	Description
12/box 85783ST	17 gauge. Standard coaxial. 0.3 mm aspiration port. Outer diameter 1.5 mm.
85784ST	17 gauge. Standard coaxial. 0.3 mm aspiration port. Outer diameter 1.5 mm. Straight shaft with 45° angled tip.
85786ST	16 gauge. Standard coaxial. 0.3 mm aspiration port. Outer diameter 1.65 mm. Straight shaft with 90° angled tip. Ideal for 12 o'clock cortical material removal.
85787ST	21 gauge. Standard biaxial. Smooth tip. 2 lateral 0.6 mm irrigation ports. Outer diameter 0.8 mm.
85794ST	1.8 mm coaxial MICS. 0.3 mm aspiration port with silicone irrigation sleeve.
85795ST	1.8 mm coaxial MICS. 0.3 mm aspiration port with silicone irrigation sleeve.
85902ST19	19 gauge. 1.8 mm biaxial MICS. 2 lateral 0.6 mm irrigation ports. Smooth tip. Outer diameter 1.07 mm.
85903ST19	19 gauge. 1.8 mm biaxial MICS. Central irrigation port. Outer diameter 1.07 mm. Sharpened inner side with atraumatic tip.
85780S	21 gauge. Standard biaxial. 0.3 mm aspiration port. Outer diameter 0.8 mm. Roughened tip.
85901S19	19 gauge. 1.8 mm coaxial MICS. 0.3 mm aspiration port. Outer diameter 1.07 mm. Roughened tip.





# Corneal Storage



## Corneal Storage

12

STERILE



Optisol-GS

Independent Corneal  
Viewing Chamber

Catalog #

Description

12/box  
50006-OPT

4°C corneal preservation medium. Permits up to 14 days preservation. Glass vials of 20mL each. Buffered tissue culture media. Medium is serum free and does not contain any bovine originating components.

IVC-12

Designed especially to reduce the possibility of tissue damage and contamination, and facilitate slit lamp and specular microscopic examination of the epithelium and endothelium. Clear viewing chamber permits easy anterior and posterior visualization of the whole cornea. Large anterior viewing window design permits wider slit lamp swing. Designed and engineered in conjunction with the Eye Bank Association of America (EBAA) and Tissue Banks International (TBI). Manufactured of durable, 100% biocompatible PMMA, to minimize potential breakage.

# Trephines



## Trephines



STERILE



### FOR IOWA P.K. PRESS PUNCH

	Catalog #	Dimensions (mm)
	1/box	Diameter
Tanne Universal Disposable Trephine Blade	E3050 6.00	6.00
Tanne Universal Disposable Trephine Blade	E3050 6.25	6.25
Tanne Universal Disposable Trephine Blade	E3050 6.50	6.50
Tanne Universal Disposable Trephine Blade	E3050 6.75	6.75
Tanne Universal Disposable Trephine Blade	E3050 7.00	7.00
Tanne Universal Disposable Trephine Blade	E3050 7.25	7.25
Tanne Universal Disposable Trephine Blade	E3050 7.50	7.50
Tanne Universal Disposable Trephine Blade	E3050 7.75	7.75
Tanne Universal Disposable Trephine Blade	E3050 8.00	8.00
Tanne Universal Disposable Trephine Blade	E3050 8.25	8.25
Tanne Universal Disposable Trephine Blade	E3050 8.50	8.50
Tanne Universal Disposable Trephine Blade	E3050 8.75	8.75
Tanne Universal Disposable Trephine Blade	E3050 9.00	9.00
Tanne Universal Disposable Trephine Blade	E3050 9.25	9.25
Tanne Universal Disposable Trephine Blade	E3050 9.50	9.50

### FOR STORZ® UNIVERSAL TREPHINE HANDLE


Trephine Blade	E3096 6.00	6.00
Trephine Blade	E3096 6.50	6.50
Trephine Blade	E3096 6.75	6.75
Trephine Blade	E3096 7.00	7.00
Trephine Blade	E3096 7.25	7.25
Trephine Blade	E3096 7.50	7.50
Trephine Blade	E3096 7.75	7.75
Trephine Blade	E3096 8.00	8.00
Trephine Blade	E3096 8.25	8.25
Trephine Blade	E3096 8.50	8.50
Trephine Blade	E3096 8.75	8.75
Trephine Blade	E3096 9.00	9.00
Trephine Blade	E3096 9.25	9.25
Trephine Blade	E3096 9.50	9.50

### FOR FREEHAND USE

Trephine Blade Long	E3096L 7.00	7.00
Trephine Blade Long	E3096L 7.25	7.25
Trephine Blade Long	E3096L 7.50	7.50
Trephine Blade Long	E3096L 7.75	7.75
Trephine Blade Long	E3096L 8.00	8.00
Trephine Blade Long	E3096L 8.25	8.25
Trephine Blade Long	E3096L 8.50	8.50
Trephine Blade Long	E3096L 8.75	8.75
Trephine Blade Long	E3096L 9.00	9.00
Trephine Blade Long	E3096L 9.25	9.25
Trephine Blade Long	E3096L 9.50	9.50



# 1.8 mm MICS™ Procedure



\* RIS – Recommended Incision Size; \* OD – Outer Diameter;  
\* ID – Inner Diameter

## VFM Phaco Packs\*



STERILE




StableChamber™ Phaco Pack  
with MICS™ Needle

StableChamber™ Phaco Pack  
with 30° Angled MICS™  
Needle

Premium Phaco Pack with  
MICS™ Needle

Premium Phaco Pack with 30°  
Angled MICS™ Needle

Catalog #	Description
6/box BL5113	1 VFM cassette assembly and I/A set with StableChamber™ tubing, 1 auxillary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 plug, 1 test chamber, 1 MICS™ needle and wrench
BL5113A	1 VFM cassette assembly and I/A set with StableChamber™ tubing, 1 auxillary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 plug, 1 test chamber, 1 angled MICS™ needle and wrench
BL5114	1 VFM cassette assembly and I/A set, 1 auxillary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 plug, 1 test chamber, 1 MICS™ needle and wrench
BL5114A	1 VFM cassette assembly and I/A set, 1 auxillary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 plug, 1 test chamber, 1 angled MICS™ needle and wrench



\* RIS – Recommended Incision Size; \* OD – Outer Diameter;  
\* ID – Inner Diameter

## AFM Phaco Packs\*



STERILE



Micro Incision AFM Phaco  
Pack with MICS™ Needle

Micro Incision AFM Phaco  
Pack with 30° Angled MICS™  
Needle

Catalog #	Description
6/box BL5123	1 AFM cassette assembly and I/A set, 1 auxillary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 test chamber, 1 MICS™ needle and wrench
BL5123A	1 AFM cassette assembly and I/A set, 1 auxillary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 test chamber, 1 angled MICS™ needle and wrench



# 2.2 mm Small Incision Procedure



## VFM Phaco Packs\*



STERILE




Premium Vacuum Phaco Pack

Basic Vacuum Phaco Pack

Optimized Stability Vacuum Phaco Pack

Catalog #	Description
6/box BL5110	1 VFM cassette assembly and I/A set, 1 auxiliary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 plug, 1 test chamber, 1 needle wrench
BL5111	1 VFM cassette assembly and I/A set, 1 irrigation sleeve, 1 plug, 1 test chamber, 1 needle wrench
BL5112	1 VFM cassette assembly and I/A set with StableChamber™ tubing, 1 auxiliary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 plug, 1 test chamber, 1 needle wrench



## AFM Phaco Packs\*



STERILE



Premium AFM Phaco Pack

Basic AFM Phaco Pack

Catalog #	Description
6/box BL5120	1 AFM cassette assembly and I/A set, 1 auxiliary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 test chamber, 1 needle wrench
BL5121	1 AFM cassette assembly and I/A set, 1 irrigation sleeve, 1 plug, 1 test chamber, 1 needle wrench



## Phaco Needles



STERILE



**MicroFlow® 2.2**  
MicroFlow® 2.2 30° Needle

**Vortex 20/21 gauge**  
Vortex 20/21 ga 30° Needle

Vortex 20/21 ga 30° Angled Needle

Catalog #	Dimensions (mm)		
6/box DP8730S	RIS*	OD*	ID*
	2.2	0.93	0.72-0.50
DP8930S	2.2	0.94 flute 0.84 shaft	0.71 flute 0.58 shaft
DP8930AS	2.2	0.94 flute 0.84 shaft	0.71 flute 0.58 shaft

\* RIS – Recommended Incision Size; \* OD – Outer Diameter; \* ID – Inner Diameter

## Required Needle Sleeve



STERILE



Required Needle Sleeve

Catalog #	Description
24/box DP5531	12 pouches, 2 per pouch

\* Also compatible with the Stellaris PC



# 2.8 mm Standard Procedure



## VFM Phaco Packs\*



STERILE




Premium Vacuum Phaco Pack

Basic Vacuum Phaco Pack

Optimized Stability Vacuum Phaco Pack

Catalog #	Description
6/box BL5110	1 VFM cassette assembly and I/A set, 1 auxiliary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 plug, 1 test chamber, 1 needle wrench
BL5111	1 VFM cassette assembly and I/A set, 1 irrigation sleeve, 1 plug, 1 test chamber, 1 needle wrench
BL5112	1 VFM cassette assembly and I/A set with StableChamber™ tubing, 1 auxiliary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 plug, 1 test chamber, 1 needle wrench



## AFM Phaco Packs\*



STERILE



Premium AFM Phaco Pack

Basic AFM Phaco Pack

Catalog #	Description
6/box BL5120	1 AFM cassette assembly and I/A set, 1 auxiliary drape, 1 screen drape, 1 remote control drape, 2 irrigation sleeves, 1 test chamber, 1 needle wrench
BL5121	1 AFM cassette assembly and I/A set, 1 irrigation sleeve, 1 plug, 1 test chamber, 1 needle wrench



## Phaco Needles



STERILE



### 19 Gauge Standard

19 ga Standard 30° Needle

19 ga Standard 45° Needle

### MicroFlow®

MicroFlow® 30° Needle

MicroFlow® 45° Needle

### MicroFlow® Plus™

MicroFlow® Plus™ 30° Needle

MicroFlow® Plus™ 45° Needle

### Vortex 19/20 Gauge

Vortex 19/20 Ga 30° Needle

Vortex 19/20 Ga 30° Angled Needle

### Thin Tip 20 Gauge

Thin Tip 20 ga 30° Needle

Thin Tip 20 ga 30° Angled Needle

Thin Tip 20 ga 45° Needle

Catalog #	Dimensions (mm)		
6/box	RIS*	OD*	ID*
DP8130S	3.0-3.2	1.7	0.9
DP8145S	3.0-3.2	1.7	0.9
DP8230S	2.65-2.75	1.07	0.90-0.51
DP8245S	2.65-2.75	1.07	0.90-0.51
DP8330S	2.75-2.85	1.17	0.98-0.61
DP8345S	2.75-2.85	1.17	0.98-0.61
DP8830S	2.65-2.85	1.16 flute 0.90 shaft	0.91 flute 0.67 shaft
DP8830AS	2.65-2.85	1.16 flute 0.90 shaft	0.91 flute 0.67 shaft
DP8430S	2.6-2.75	0.90	0.67
DP8430AS	2.6-2.75	0.90	0.67
DP8445S	2.6-2.75	0.90	0.67

\* RIS - Recommended Incision Size; \* OD - Outer Diameter; \* ID - Inner Diameter

## Required Needle Sleeve



STERILE



Required Needle Sleeve

Catalog #	Description
24/box DP5531	12 pouches, 2 per pouch

# Other Cataract Products



## Anterior Vitrectomy Packs



STERILE



20 Gauge Pneumatic Anterior Cutter with Irrigation

Catalog #

6/box

DP4400-6

20 Gauge Pneumatic Anterior Cutter without Irrigation

MVS3201-6



## Bipolar Accessories



STERILE



Stellaris® Bipolar Cord, Lemo Connector

Catalog #

12/box

CX9430

Straight Bipolar Forceps  
0.5 mm Tip with Lemo Connector

D8200



Bipolar Scleral Pencil, 45° Beveled Tip with Lemo Connector

E7928



## DigiFlow™



STERILE



DigiFlow™ Tubing

Catalog #

6/box

D4600A

Description

For use with the DigiFlow software upgrade which enables control of automatic pressurization of the BSS bottle, minimizing IOP fluctuations during cataract surgery.



## Maintenance



Universal Maintenance Kit

Catalog #

1/box

CX7120



## Tray



Stellaris® Ultrasonic Handpiece Tray

Catalog #

1/box

BL3379



# 23 Gauge Procedure

## Procedure Packs



STERILE



**STELLARIS™ PC**  
23 Gauge Posterior  
Procedural Pack

Catalog #

6/box  
BL5223

Description

23 Gauge Combined  
Procedural Pack

BL5323

23 Gauge Posterior  
Procedural Pack with  
Integrated Vented AFI

BL5423

23 Gauge Combined  
Procedural Pack with  
Integrated Vented AFI

BL5523



## MILLENNIUM® MVE

Premium TSV23™ MVE  
Vitrectomy Pack - 2500 CPM

USSCX6923A

Preassembled Pack Includes:

- MVE 23G High Speed Cutter
- Venturi Cassette With Stopper & Cap
- Irrigation Administration Tube Set
- Filtered Air Exchange Line
- Utility Line and Male-Male Connector
- 10cc Syringe
- Mayo tray, Remote Control and Screen Drapes
- (2) – 3-Way Stopcocks
- Green Pneumatic Actuation Tubing
- (3) – Trocar/Cannula Assemblies
- 23G 5mm Infusion Line
- (3) – 23G Cannula Plugs
- TSV23™ Light Pipe

# 23 Gauge Procedure



6

2

STERILE



6

STERILE



2

STERILE



8

STERILE



6

STERILE



## Instrumentation

### FORCEPS

23 Gauge Single-Use Endgripping Forceps

23 Gauge Single-Use Asymmetrical Peeling Forceps

### KNIVES

23 Gauge MVR Knife

### LIGHT PIPE

TSV23™ Light Pipe

### BACKFLUSH HANDPIECE

23 Gauge Backflush Handpiece with Soft Tip Cannula

### MEMBRANE PICK

23 Gauge Membrane Pick - 130°

### CANNULAS

23 Gauge Single Trocar/Cannula  
23 Gauge Soft Tip Cannula

Catalog #

Description

6/box  
65023S

2/box  
65024S

Ultra-thin endgripping tip for enhanced visibility and grasp. Ergonomic handle consistent with our popular reusable design. Balanced design for optimal feel and control.

SUI03

SUI04

Used for ILM and Epiretinal Membrane Peel. Use on macular pickers and cellophane maculopathy. Asymmetric jaw allows for excellent visualization of retinal tissue

6/box  
E7523

A=31.8 mm; B=0.65 mm; C=1.1 mm

6/box  
CX9623

2/box  
SUJ04

2 per box.

8/box  
SUK01

8 per box.

6/box  
CX9123

10/box

6 per box.

BL7823

Aspirating cannula with soft silicone tip. 10 per box.



# 25 Gauge Procedure

## Procedure Packs



STERILE



### STELLARIS™ PC

25 Gauge Posterior  
Procedural Pack

Catalog #

6/box  
BL5225

Description

25 Gauge Combined  
Procedural Pack

BL5325

25 Gauge Posterior  
Procedural Pack with  
Integrated Vented AFI

BL5425

25 Gauge Combined  
Procedural Pack with  
Integrated Vented AFI

BL5525

### MILLENNIUM® MVE

Premium TSV25™ MVE  
Vitrectomy Pack - 2500 CPM

USSCX6925A

Preassembled Pack Includes:

- MVE 25G High Speed Cutter
- Venturi Cassette With Stopper & Cap
- Irrigation Administration Tube Set
- Filtered Air Exchange Line
- Utility Line and Male-Male Connector
- 10cc Syringe
- Mayo tray, Remote Control and Screen Drapes
- (2) – 3-Way Stopcocks
- Green Pneumatic Actuation Tubing
- (3) – Trocar/Cannula Assemblies
- 25G 5mm Infusion Line
- (3) – 25G Cannula Plugs
- TSV25™ Light Pipe

### MILLENNIUM®

TSV25™ Electric  
Vitrectomy Pack - 1500 CPM

USSCX4925A

Pack Includes:

- 25G High Speed Electric Cutter
- Venturi Cassette With Stopper & Cap
- Irrigation Administration Tube Set
- Filtered Air Exchange Line
- Utility Line and Male-Male Connector
- 10cc Syringe
- Mayo Tray, Remote Control and Screen Drapes
- (2) – 3-Way Stopcocks
- (3) – Trocar/Cannula Assemblies
- 4mm Infusion Cannula
- (3) – 25G Cannula Plugs
- TSV25™ Light Pipe



# 25 Gauge Procedure



6

2

STERILE



6

STERILE



2

STERILE



8

STERILE



10

STERILE



## Instrumentation

### FORCEPS

25 Gauge Single-Use Endgripping Forceps

25 Gauge Single-Use Asymmetrical Peeling Forceps

### KNIVES

25 gauge MVR Knife

TSV25™ ACCESSORIES  
25 Gauge Tapered Endolase Probe

TSV25™ Light Pipe

TSV25™ Aspirating Pick

TSV25 Extendable Pick (Blunt Tip)

### BACKFLUSH HANDPIECE

25 Gauge Backflush Handpiece with Soft Tip Cannula

### MEMBRANE PICK

25 Gauge Membrane Pick - 130°

### CANNULAS

25 Gauge Soft Tip Cannula

Catalog #

Description

6/box  
65025S

2/box  
65026S

Ultra-thin endgripping tip for enhanced visibility and grasp. Ergonomic handle consistent with our popular reusable design. Balanced design for optimal feel and control.

SUI05

SUI06

Used for ILM and Epiretinal Membrane Peel. Use on macular pucker and cellophane maculopathy. Asymmetric jaw allows for excellent visualization of retinal tissue

6/box  
E7525

A=31.8 mm; B=0.5 mm; C=1.0 mm

6/box  
CX5985

CX9625

CX9627

CX9678

2/box  
SUJ05

2 per box.

8/box  
SUK02

8 per box.

10/box  
BL7725

Aspirating cannula with soft silicone tip.



# 20 Gauge Procedure

## Procedure Packs



STERILE



**STELLARIS™ PC**  
20 Gauge Posterior  
Procedural Pack

Catalog #

6/box  
BL5220

Description

20 Gauge Combined  
Procedural Pack

BL5320

20 Gauge Posterior  
Procedural Pack with  
Integrated Vented AFI

BL5420

20 Gauge Combined  
Procedural Pack with  
Integrated Vented AFI

BL5520



## MILLENNIUM® MVE

Premium 20 Gauge MVE  
Vitrectomy Pack – 2500 CPM

USSCX6920A

Preassembled Pack Includes:

- MVE 20G High Speed Cutter
- Venturi Cassette With Stopper & Cap
- Irrigation Administration Tube Set
- Filtered Air Exchange Line
- Utility Line and Male-Male Connector
- 10cc Syringe
- Mayo Tray, Remote Control and Screen Drapes
- (2) – 3-Way Stopcocks
- Green Pneumatic Actuation Tubing
- Scleral Plugs (2 each in 19G & 2 20G)
- 20G 4mm Infusion Cannula
- 20G MVR Blade (Straight)
- 20G Light Pipe



# 20 Gauge Procedure



## Procedure Packs



STERILE



### MILLENNIUM®

20 Gauge Pneumatic  
Vitrectomy Pack - 750 CPM

20 Gauge Electric Vitrectomy  
Pack - 1500 CPM

Catalog #

Description

6/box

USSCX4820A

Pack Includes:

- 20G Pneumatic Cutter
- Venturi Cassette With Stopper & Cap
- Irrigation Administration Tube Set
- Filtered Air Exchange Line
- Utility Line and Male-Male Connector
- 10cc Syringe
- Mayo Tray, Remote Control and Screen Drapes
- (2) - 3-Way Stopcocks
- Scleral Plugs (2 each in 19G & 2 20G)
- 4mm Infusion Cannula
- 20G MVR Blade (Straight)
- 20G Light Pipe

USSCX4920A

Pack Includes:

- 20G High Speed Electric Cutter
- Venturi Cassette With Stopper & Cap
- Irrigation Administration Tube Set
- Filtered Air Exchange Line
- Utility Line and Male-Male Connector
- 10cc Syringe
- Mayo Tray, Remote Control and Screen Drapes
- (2) - 3-Way Stopcocks
- Scleral Plugs (2 each in 19G & 2 20G)
- 4mm Infusion Cannula
- 20G MVR Blade (Straight)
- 20G Light Pipe



## Instrumentation/Accessories



STERILE



### KNIVES

20 Gauge MVR Knife

### LIGHT PIPES

20 Gauge Cannonball Wide  
Angle Light Pipe

20 Gauge Light Pipe

### NEEDLE/CUTTER

20 Gauge Fragmentation  
Needle 15°

20 Gauge High Speed  
Anterior Vitrectomy Cutter

Catalog #

Description

6/box

E7520

A=31.8 mm; B=1.2 mm;  
C=2.4 mm

DP9604

MVS1011C

BL5275












BL5612





# Other Retina

## Non-Gauge Specific Accessories

		Catalog #	Description
	 <b>STERILE</b> 	6/box D4600A	Air Exchange Line
		CX5710	<b>VISCOUS FLUID</b> Viscous Fluid Pack For Silicone Aspiration
		CX5720	Viscous Fluid Pack For Injection of ES 5000 S Silicone Oil Syringe
	 <b>STERILE</b> 	1/box ES 5000 S	10ml syringe. To be used with CX5720 Viscous Fluid Pack.
	 <b>STERILE</b> 	10/box TR9060 61	<b>INTRAOCULAR GAS</b> Universal Gas Kit



## 2 Pin Connection

10

### 18 GAUGE

18 Gauge Straight

Catalog #

Description

E7918

10 per box.



STERILE

18 Gauge Curved

E7918 C

10 per box.



### 20 GAUGE

20/23 Gauge Tapered

E7923

10 per box.



12

### 23 GAUGE

23 Gauge Tapered

E7931

Designed for retinal endodiatheirmy. 12 per box.

STERILE



### 25 GAUGE

25 Gauge Tapered

E7924

Designed for retinal endodiatheirmy. 12 per box.



25 Gauge Blunt

E7925

Designed for retinal endodiatheirmy. 12 per box.



10

### CORD

Disposable Bipolar Cord

S2070 C

10 per box.

STERILE



1

### ADAPTER

Reusable Banana Plug Adapter

CX9404

Reusable.



## Lemo Connection

12

### 23 GAUGE

23 Gauge Tapered

Catalog #

Description

E7930

Designed for retinal endodiatheirmy. 12 per box.

STERILE





# Crystalens® Accommodating IOL



## Crystalens Accommodating IOL



STERILE



	Catalog #	Available Powers	Shape	Material	Diameter (mm)	Rec. Starting A Constant	Rec. Starting ACD (mm)	Overall Length (mm)
Crystalens AO™	AT50AO	17.0-33.0 D in 0.50-D steps	Biconvex	Biosil®	5.0	119.1	5.61	11.5
	AT52AO	18.0-22.0 D in 0.25-D steps	Biconvex	Biosil®	5.0	119.1	5.61	12.0
Crystalens HD™	HD500	10.0-16.5 D in 0.50-D steps	Biconvex	Biosil®	5.0	118.8	5.43	11.5
	HD520	17.0-33.0 D in 0.50-D steps	Biconvex	Biosil®	5.0	118.8	5.43	12.0
Crystalens Five-O™	AT50SE	18.0-22.00 D in 0.25-D steps	Biconvex	Biosil®	5.0	119.0	5.55	11.5
	AT52SE	27.0 to 33.0 D in 0.50-D steps	Biconvex	Biosil®	5.0	119.0	5.55	12.0



## Crystalens Delivery System




STERILE



Crystalsert™

Catalog #	Description
CI-28	Delivery system for Crystalens AO™, Crystalens HD™ and Crystalens Five-O™


# Akreos® Monofocal IOL



## Akreos Monofocal IOL\*



STERILE



		Available Powers	Shape	Material	Optic Size (mm)	Overall Length (mm)	Applanation A Constant ACD Constant Surgeon Factor	IOL Master** A Constant ACD Constant Surgeon Factor
Akreos MICS™	MI60L	0.0-9.0 D in 1.0-D steps	Asymmetrical Biconvex	Hydrophilic Acrylic (26% Water Content)	5.6 22.5-30.0D	10.5 22.5-30.0D	118.4	118.9
		10.0-30.0 D in 0.5-D steps	4 Haptics at 8° Angle		6.0 15.5-22.0D	10.7 15.5-22.0D	5.2 mm 1.45 mm	5.49 mm 1.73 mm
					6.2 0.0-15.0D	11.0 0.0-15.0D		
Akreos AO™	AO60	0.0-9.0 D in 1.0-D steps	Asymmetrical Biconvex	Hydrophilic Acrylic (26% Water Content)	5.6 22.5-30.0D	10.5 22.5-30.0D	118.0	118.5
		10.0-30.0 D in 0.5-D steps	4 Haptics at 0° Angle		6.0 15.5-22.0D	10.7 15.5-22.0D	4.96 mm 1.22 mm	5.18 mm 1.40 mm
					6.2 0.0-15.0D	11.0 0.0-15.0D		

\*UV Absorbing; 360° Posterior Square Edge; Refractive Index: 1.458

\*\* Updated on 10/2009. Subject to future updates.



## Akreos Lens Inserters



STERILE




		For Lens Model	Recommended Incision Size*** (mm)	Type of Action	Material
Viscoject 1.8	LP604350	MI60L	1.8-2.4	Syringe	Plastic
		AO60	2.8-3.0	Syringe	Polypropylene
AI-28	AI-28	AO60	2.8-3.0	Syringe	Polypropylene
Viscoject 2.2	LP604340	AO60	2.2-2.8	Syringe	Plastic

\*\*\*Incision sizes are estimates. Actual sizes may vary due to incision construction, location and any stretching that may have occurred.



# SofPort® IOL



## SofPort IOL\*



STERILE



	Catalog #	Available Powers	Shape	Material	Optic Size (mm)	Overall Length (mm)	Applanation A Constant ACD Constant Surgeon Factor	IOL Master A Constant ACD Constant Surgeon Factor
SofPort AO	LI61AO	0.0-4.0 D in 1.0-D steps	Asymmetrical Biconvex	Advanced Generation Silicone	6.0	13.0	118.0	118.7
		5.0-30.0 D in 0.5-D steps	Modified C Haptics at 5° Angle				5.0 mm 1.22 mm	5.40 mm 1.62 mm
SofPort	LI61SE	0.0-4.0 D in 1.0-D steps	Equiconvex	Advanced Generation Silicone	6.0	13.0	118.0	118.7
		5.0-30.0 D in 0.5-D steps	Modified C Haptics at 5° Angle				5.0 mm 1.22 mm	5.40 mm 1.62 mm

\*360° Posterior Square Edge; Blue PMMA Haptics; Refractive Index: 1.43



## SofPort System Inserters



STERILE



	Catalog #	For Lens Model	Recommended Incision Size** (mm)	Type of Action	Material
EZ-28	EZ-28	LI61AO LI61SE	2.8-3.0	Syringe	Polypropylene

\*\*Incision sizes are estimates. Actual sizes may vary due to incision construction, location and any stretching that may have occurred.

# PMMA Lenses



## Single-Piece PMMA Lenses



STERILE



	Catalog #	Available Powers	Shape	Haptics	Optic Size (mm)	Overall Length (mm)	Applanation A Constant Surgeon Factor	Other Features
EZE-50	EZE-50	10.0-30.0 D in 0.5-D steps	Equiconvex	Modified C Double Flat Step-Vaulted 3° Angle	5.0	12.75	118.1 5.02 mm 1.28 mm	EZE-FIT® 0
EZE-55	EZE-55	10.0-30.0 D in 0.5-D steps	Equiconvex	Modified C Double Flat Step-Vaulted 3° Angle	5.5	12.75	118.1 5.02 mm 1.28 mm	EZE-FIT® 0
P359UV	P359UV	0.0-45.0 D in 1.0-D steps 10.0-30.0 D in 0.5-D steps	Equiconvex	Modified C Double Flat Step-Vaulted 2.4° Angle	5.5	12.25	118.0 4.96 mm 1.22 mm	EZE-FIT® 2 Extended diopter range available
EZE-60	EZE-60	10.0-30.0 D in 0.5-D steps	Equiconvex	Modified C Double Flat Step-Vaulted 3° Angle	6.0	12.75	118.1 5.02 mm 1.28 mm	EZE-FIT® 0
P574UV	P574UV	-18-9.0 D in 1.0-D steps	Equiconvex	Modified C Double Flat Step-Vaulted 3° Angle	6.0	13.25	118.1 5.02 mm 1.28 mm	EZE-FIT® 2 Extended diopter range available
EZE-65	EZE-65	10.0-30.0 D in 0.5-D steps	Equiconvex	Modified C Double Flat Step-Vaulted 3° Angle	6.5	13.25	118.1 5.02 mm 1.28 mm	EZE-FIT® 0
P366UV	P366UV	10.0-30.0 D in 0.5-D steps	Equiconvex	Modified C Step-Vaulted 6.6° Angle	6.5	13.4 12.0 eyelet to eyelet	118.5 5.26 mm 1.51 mm	Haptic Eyelets 0.43 mm
EZE-70	EZE-70	10.0-30.0 D in 0.5-D steps	Equiconvex	Modified C Double Flat Step-Vaulted 3° Angle	7.0	13.5	118.1 5.02 mm 1.28 mm	EZE-FIT® 0
S122UV	S122UV	5.0-30.0 D in 0.5-D steps	Equiconvex	Modified S Step-Vaulted 4.4° Angle	6.0	12.5	115.8 3.68 mm -0.02 mm	Anterior Chamber 4 Point Fixation <u>Sizing Guidance</u> White to White Range 10.0 mm - 11.49 mm (Corneal Diameter)
L122UV	L122UV	5.0-30.0 D in 0.5-D steps	Equiconvex	Modified S Step-Vaulted 3.7° Angle	6.0	13.75	115.8 3.68 mm -0.02 mm	Anterior Chamber 4 Point Fixation <u>Sizing Guidance</u> White to White Range 11.5 mm - 12.25 mm (Corneal Diameter)

The EZE-FIT IOL series of lenses provide a choice of haptic flexibility. EZE-FIT 2 = standard flexibility haptics; EZE-FIT 0 = more flexible haptics.



# Viscoelastics



## Sodium Hyaluronate



STERILE



Amvisc® Plus 0.5mL

Amvisc® Plus 0.8mL

Amvisc® 0.5mL

Amvisc® 0.8mL

Catalog #	Description
1/box 60051	Enhanced high viscosity pure sodium hyaluronate. 1.6% concentration. 0.5mL preloaded sterile glass syringe with Luer-Lok and 27 gauge disposable cannula. Minimum order of 5.
60081	Enhanced high viscosity pure sodium hyaluronate. 1.6% concentration. 0.8mL preloaded sterile glass syringe with Luer-Lok and 27 gauge disposable cannula. Minimum order of 5.
59051	High viscosity pure sodium hyaluronate. 1.2% concentration. 0.5mL preloaded sterile glass syringe with Luer-Lok and 27 gauge disposable cannula. Minimum order of 5.
59081	High viscosity pure sodium hyaluronate. 1.2% concentration. 0.8mL preloaded sterile glass syringe with Luer-Lok and 27 gauge disposable cannula. Minimum order of 5.

## 2% Hydroxypropylmethylcellulose



STERILE



Ocucoat® Viscoadherent



STERILE



Ocucoat® Viscoadherent


Catalog #	Description
1/box CC050	2% HPMC viscoadherent solution. 1.0mL preloaded sterile glass syringe with Luer-Lok and 25 gauge disposable cannula. 1 per box.
6/box CC065	2% HPMC viscoadherent solution. 1.0mL preloaded sterile glass syringe with Luer-Lok and 25 gauge disposable cannula. 6 per box.



# Eye Shields




## Eye Shields



12

STERILE



Universal Eye Shield, Adult

Universal Eye Shield, Child

Universal Goffman Blue Eye Garter Shield Cover and Blue Eye Shield

Catalog #	Dimensions (mm)	Description
12/box E5689 A	60 x 75	Clear plastic. Accommodates either eye.
E5689 B	53.8 x 69	Clear plastic. Accommodates either eye.
E5690 U	60.3 x 76.2	Blue plastic shield for reduced glare.
50/box E5692 D	60.3 x 76.2	Accommodates either eye. Slotted malleable aluminum surface. Correct size holes produce ideal pin hole effect.
E5691 A50	60.3 x 76.2	Assorted colors. Includes (1) E5692 D aluminum Fox Eye Shield and (1) Eye Garter Shield Cover.
E5691 W50	60.3 x 76.2	White only. Includes (1) E5692 D aluminum Fox Eye Shield and (1) Eye Garter Shield Cover.

50

STERILE




Fox Eye Shield

Combination Eye Garter Shield Cover and Fox Shield, Assorted Colors

Combination Eye Garter Shield Cover and Fox Shield, White Only

## Eye Shields Accessories



50

STERILE



Eye Garter Shield Cover, Assorted Colors

Eye Garter Shield Cover, White Only

100

STERILE




STORZ® Ophthalmics Easy Shield Adhesive Bandages

Catalog #	Dimensions (mm)	Description
50/box E5693		Assorted colors. Used with E5692 D and E5697 A eye shields to cushion the eye. Eliminates need for taping. No sticky residue.
E5693 W		White only. Used with E5692 D and E5697 A eye shields to cushion the eye. Eliminates need for taping. No sticky residue.
100/box E5691 E		24 packages of 4 each. For E5691 series eye shield sets. Comfortable 3M hypoallergenic medical foam adhesive.



# Corneal Shields & Other




## Corneal Shields



Crouch Corneal Protector,  
Adult

Catalog #	Dimensions (mm)	Description
10/box E5699	13.5 x 25.8	For oculoplastic procedures for corneal and scleral protection. Dimensions enable protector to rest on the sclera and not the cornea. Easily removed with E1500 Bishop-Harmon Forceps.
E5699 P	17	For pediatric oculoplastic procedures for corneal and scleral protection. Dimensions enable protector to rest on the sclera and not the cornea. Easily removed with E1500 Bishop-Harmon Forceps.
6/box G11X-S	16	Dissolves in approximately 12 hours.
G145X-S	14.5	Dissolves in approximately 12 hours.



Pediatric Corneal Shield



SurgiLens® Corneal Shield  
SurgiLens® Corneal Shield



## Other Accessories



Fluid Filter

Catalog #	Description
50/box E4429 22	Flow rate 45cc/minute. Pore size 0.2 micron. Bi-directional membrane support rated at 5.3kg/cm <sup>2</sup> (75 psi) via prefilled syringe and aspiration from vials.
E4429 45	Flow rate 80cc/minute. Pore size 0.45 micron. Bi-directional membrane support rated at 5.3kg/cm <sup>2</sup> (75 psi) via prefilled syringe and aspiration from vials.
1/box E5675	Non-reflective black anodized aluminum. Selection of five pin hole sizes (0.8mm to 1.5mm). Overall length 102mm, 4.0 inches.



Goffman Occluder



Balanced Salt Solution



Balanced Salt Solution

BS018	18mL.
BS050	500mL.

