

Product	Attributes	ISO 11607 Part 1	EN868 Part 2
Safeguard Regular Green/blue interleaved Green only Blue only	An ideal general usage sterilisation wrap <ul style="list-style-type: none"> • 57 gsm • Alcohol repellent • Medium weight and strength • High burst strength 	☑	☑
Safeguard Extra Green/blue interleaved	An excellent barrier to aqueous borne bacteria <ul style="list-style-type: none"> • 67 gsm • Water and alcohol repellent • Good strength • Tear resistant 	☑	☑
Safeguard 50 Green/blue interleaved	Combines a strong blue outer wrap which has a high filtration medium with an alcohol repellent green inner non-woven <ul style="list-style-type: none"> • 57 gsm inner green wrap • 50 gsm outer blue wrap • Toughness combined with economy • High burst strength 	☑	☑
Safeguard 100 Green/blue interleaved	A totally superior tray wrap. It has physically strong blue and green wraps that provide a high filtration medium. Where consistent tearing is a problem, can help reduce expensive recycling <ul style="list-style-type: none"> • 40 gsm inner green wrap • 50 gsm outer blue wrap • Lint resistant • High filtration • Toughness combined with durability 	☑	☑
RHINO WRAP Blue outer and green inner wraps, bonded together to create a really tough single wrap	Double protection with a single action, saving time but allowing secure, safe wrapping of trays <ul style="list-style-type: none"> • 80 gsm • Toughness combined with durability • Reduced time - increased efficiency 	☑	☑
Crepe paper White, blue or green	Sterilisation wrapping <ul style="list-style-type: none"> • 60 gsm • Suitable for steam, ETO and Gamma sterilisation 	☑	☑
Plain paper White, blue or green	Sterilisation wrapping <ul style="list-style-type: none"> • 60 gsm • Suitable for steam, ETO and Gamma sterilisation 	☑	☑
Tray Liner White only	Highly absorbent, nonwoven pad <ul style="list-style-type: none"> • 80 gsm • Absorption capacity 225% • Maintains strength & uniformity • Suitable for steam sterilisation 	☑	☑
Transport Wrap Purple	Protective wrap - not intended to form part of the sterilisation process <ul style="list-style-type: none"> • 60 gsm • Provides extra protection during the transportation of wrapped trays 	N/A	N/A

Uncontrolled copy – Accurate at time of issue

Doc No	QA Date
QA 417/4	03/10/11



BDF – Grangestone - GIRVAN – Ayrshire – KA26 9PF
Guardian line: 01465 716716 Fax: 01465 714438
www.bdf.ltd.uk sales@bdf.ltd.uk

BDF trades as part of Sunlight Clinical Solutions, a division of The Sunlight Service Group Ltd