

# Disposable Polypropylene Hollowware





The Whelehan Group offers a comprehensive range of disposable polypropylene hollowware in natural, pink, blue or green.

Our wide range of products comprises eight items in the Gallipot and Bowl range, three sizes of Receiver and three Jugs. We also offer eight Tray options and a 20cm Quiver. These are compatible with Autoclave, Gamma or Ethylene Oxide sterilisation. Products are all available sterile or non-sterile, single or double wrapped, or pre assembled into bespoke packs. Some sterile items are available in a peel pouch presentation.

All products are latex free.

Tel: 353 | 8068600 Fax: 353 | 836227 | Email: sales@tpwhelehan.ie www.whelehangroup.com

3 Northern Cross Business Park Finglas, Dublin 11 Ireland



Receivers



Bowls



Quivers

## **Individual Products**



Trays



AVM 4022 Maternity (Placenta) Bowl Set



Jugs



AVM 4031 Pink Polyware Set



Gallipots



AVM 4042 Basic Polyware Pack

## **Pre Assembled Packs**



AVM 4023 Laparotomy Bowl Set



AVM 4026 Basic Set



AVM 4040 Theatre Bowl



AVM 4041 Bowl Set Large



AVM 4044 Supplementary Polyware Set



AVM 4043 Bowl Set Small





Products are available vacuum packed to improve nesting and reduce storage space





Easy opening peel apart packaging





Specialist Tyvek pouches allow EO Sterilisation





Bowl Liners available bulk packed or pre sterilised

### Individual Product Details

Non Sterile Code	Product Description	Case Qty	Case Dimensions (mm)
AVM0146	Kidney Dish (8'') Graduated 200×100×40mm	380	476×233×410
AVM0143	Kidney Dish (10") Graduated 250x120x50mm 66g	100	489×359×296
AVM0243	Kidney Dish (10") Graduated 250x120x50mm 46g	200	511×331×455
AVM0139	Kidney Dish (12") Graduated 300x150x60mm	150	360×360×430
AVM0328	Tray 2 compartment 183×100×25mm	500	500×325×430
AVM0314	Tray 3 compartment 183×100×25mm	500	500×325×430
AVM0316	Tray 4 compartment 183×100×25mm	500	500x325x430
AVM0327	Tray with 1 integral gallipot 196×130×28mm	200	582×223×295
AVM0315	Tray with 2 integral gallipots 196×130×28mm	200	595×223×295
AVM0326	Tray 196x130x28mm	200	582×223×295
AVM0318	Tray 4 section 270x180x42mm	52	650x286x370
AVM0312	Tray 300x250x50mm	50	643×277×325
AVM0320	Tray 4000ml 290x220x90mm	60	600×330×500
AVM0396	Gallipot 60ml Graduated	2000	500×290×330
AVM0395	Gallipot 120ml Graduated	1000	496×331×350
AVM0340	Bowl / Gallipot 250ml Graduated	750	511x331x455
AVM0325	Bowl 500ml Graduated 119mm diameter x 65mm	100	496×331×350
AVM0321	Bowl 1000ml Graduated 170mm diameter x 80mm	400	496×331×350
AVM0338	Bowl 2000ml 200mm diameter x 95mm Graduated	200	496×331×350
AVM0333	Bowl 4000ml (10") 250mm diameter x 130mm Graduated	100	476×233×410
AVM0354	BowlTheatre (14")	50	360x360x430
AVM0364	Theatre bowl liner	50	360×360×430
AVM0344	Jug / Bowl 500ml Graduated	40	582×223×295
AVM0385	Jug 500ml Graduated	200	496×331×350
AVM0427	Jug 1000ml Graduated	100	511x331x455
AVM0460	Diathermy Quiver 20cm	50	360×360×430
AVM0383	Medicine Measure 25ml Graduated	4000	489×359×296
AVM0382	Medicine Measure 30ml Graduated	2000	415×257×219
AVM0387	Medicine Measure 50ml Graduated	1000	415×257×219
AVM0384	Medicine Measure 60ml Graduated	1000	415×257×219
AVM0386	Medicine Spoon 5ml Graduated	2000	495×266×266

### AVM4022 Maternity (Placenta) Bowl Set

1 x Theatre Bowl 5 litre 1 x Kidney Dish 250mm (10") 46g

I × Bowl 2 litre
I × Bowl 500ml
I × Bowl 250ml
Wrapped in I20 × I20cm PETissue and Polythene Cover

### AVM4040 Theatre Bowl

I x Theatre Bowl 5 litre Pink Wrapped in 120 x 120cm PETissue and Polythene Cover

AVM4042 Basic Polyware Pack 2 x Kidney Dish 300mm (12'') 3 x Bowl 250ml 1 x Kidney Dish 200mm (8") I x Diathermy Quiver Wrapped in 120 x 120cm PETissue and Peel Pouch

### AVM4023 Laparotomy Bowl Set

I × Theatre Bowl 5 litre 2 × Kidney Dish 250mm (10'') 46g I x Jug I litre 3 x Bowls 250ml I × Tray 2 gallipot Double wrapped in I20 × I20cm PE Tissue and Polythene Cover

### AVM4031 Pink Polyware B Set

3 x Kidney Dish 250mm (10") 46g Pink 2 x Bowl 250ml Pink Double wrapped in 60 x 60cm PETissue and Polythene Cover

### AVM4043 Bowl Set Small

2 x Kidney Dish 250mm (10") 66g 2 x Bowl 250ml 1 x Bowl 1 litre I x Diathermy Quiver Wrapped in 90 x 90cm Green Paper and Polythene Cover

### AVM4026 Basic Set

2 x Kidney Dish 250mm (10") 46g  $3 \times Bowl 250ml$ I x Diathermy Quiver Double wrapped in 60 x 60cm PETissue and Polythene Cover

### AVM4041 Bowl Set Large

3 × Kidney Dish 250mm (10") 66g 2 × Bowl 250ml 1 x Jug 1 litre I x Diathermy Quiver Wrapped in 90 x 90cm Green Paper and Polythene Cover

### AVM4044 Supplementary Polyware Set

2 × Kidney Dish 300mm (12") 2 × Bowl 250ml Wrapped in 120 × 120cm PETissue and Peel Pouch