

ANTERIOR CHAMBER INTRAOCULAR LENSES

S122UV L122UV

S122UV

Optic: 6.0mm, Equiconvex PMMA

Haptics: Step-Vaulted
4.4° Effective Angle
4 Point Fixation

Length: 12.50mm

Diopters: +10 to +30

Lens Constants:

A-Constant: 115.8

ACD: 3.68mm

Surgeon Factor: -0.02

L122UV

Optic: 6.0mm, Equiconvex PMMA

Haptics: Step-Vaulted
3.7° Effective Angle
4 Point Fixation

Length: 13.75mm

Diopters: +5 to +30


Lens Constants:

A-Constant: 115.8


ACD: 3.68mm

Surgeon Factor: -0.02

S122UV



L122UV



NOTE: A-Constant, ACD, and Surgeon Factor are estimates only. It is recommended that each surgeon develop his or her own values.

CAUTION: Federal (U.S.A.) law restricts this device to sale by, or on order of, a physician.

Please see package insert for indications, contraindications, and precautions prior to use.

© 2000 Bausch & Lomb Surgical, Inc.

All Rights Reserved

**BAUSCH
& LOMB**
Surgical

1-800-338-2020
<http://www.blsurgical.com>

SPA-6903