



Fine Surgical Instruments



Ophthalmic





the right products
in the right hands
at the right price

To provide our customers with the service and expertise they would expect from a manufacturer of such critically important products used in the day-to-day running of their clinical department.





what makes Unisurge different?

Established in 1988, Unisurge is a trusted manufacturer and supplier of practitioner-specified procedure packs, dressing packs, medical disposable and theatre products to health care professionals. Our range of bespoke single use products is both cost effective and patient specific. Surgical teams can be confident that all the necessary components for a procedure will be supplied, to their specification, in a sterile pack, ready for immediate use.

Our promise to you:

Experienced sales team ready to take your order and advise you on delivery times and other Unisurge products – we aim to surpass all the standards of your current service.

Innovative products constantly being introduced to increase your purchasing efficiency.

A comprehensive range of products always kept in stock ready for fast-track delivery to your door.

Regular deliveries throughout the UK.

All at the most cost effective prices.



our commitment to you...

From our Newmarket headquarters, the development team keeps abreast of the latest advances in surgical practice and takes full account of the consistent rise in patient expectations.

Unisurge offers a comprehensive range of theatre quality surgical instruments. They are available as re-usable as well as single-use products. These include individual supplementary instruments, instrument sets and instruments that can be supplied as part of a sterile procedure pack which can include all the components for any surgical procedure regardless of the speciality. The contents of each of our packs can be practitioner specified to meet the requirement of all modern surgical procedures including drapes and gowns.

Linked to a significant investment in state of the art information technology, we provide a proven and trusted distribution network. Using our own fleet of vehicles, Unisurge is able to ensure that sterile packs arrive at their destination on time and in optimum condition.


**Our aim is to tailor our service to suit your specific needs
and quite simply surpass the expectations you place upon us.**

Signature Instruments...

A comprehensive range of instruments that have been manufactured and hand finished by craftsmen. There might be slight variations in the dimensions of an instrument from those given in the descriptions however, the instruments are inspected to ensure that any variations are within agreed tolerances. Unisurge operates active Quality Control / Assurance to fully comply with the regulatory requirements.

All instruments are manufactured from stainless steel and have been finished to the highest standards.

Unisurge is committed to achieving full compliance relating to quality and labour standards.





other information



Single use. Do not re-use



Lot/Batch number



Date of manufacture

REF

Replaces the use of the words such as code, reference, article number



Sterile. Sterilised by Ethylene Oxide.



Read instructions for use and/or warnings included



Shows "use-by" date for the product as YYYY-MM



Latex-free

Introducing our Ophthalmic Range



Fine Surgical Instruments

Type	Ref No	Description	Case Quantity
Clamps			
	Chalazion		
	F895009	Lambert 10mm Round Plate	12 X 20
	F895020	Lambert 15mm Round Plate	12 X 20
	F895013	Desmarres 13mm Oval Plate	12 X 20
	F895010	Desmarres 16mm Oval Plate	12 X 20
	F895011	Desmarres 18mm Oval Plate	12 X 20
	F895012	Desmarres 20mm Oval Plate	12 X 20
	F895067	Desmarres 31mm Oval Plate	12 X 20
	F895030	Greene Tarsal Cyst 13mm Round	12 X 20
Curettes			
	Spoon		
	F895093	Lambert Size A Round 1.5mm	12 X 20
	F895068	Lambert Size B Round 2mm	12 X 20
	F895002	Lambert Size C Oval 3 x 2.2mm	12 X 20
	F895003	Lambert Size D Oval 3.8 x 2.7mm	12 X 20
	F895004	Lambert Size E Oval 4.6 x 3mm	12 X 20
	F895005	Lambert Size F Oval 5.7 x 3.7mm	12 X 20
	F895066	Meibomian Cyst Double Ended 1mm & 2mm	12 X 20
Dilators			
	Lacrimal		
	F895007	Foster	12 X 20
	F895008	Nettleship	12 X 20
Handles			
	F893002	Mini Blade	12 X 20
Hooks			
	Strabismus		
	F895031	Moorfields 14cm (5.6")	12 X 20
	F895077	Kilner (Gillies) 14cm (5.6")	12 X 20
	F895073	Graefe Small 14cm (5.6")	12 X 20
	F895072	Graefe Medium 14cm (5.6")	12 X 20
	F895071	Graefe Large 14cm (5.6")	12 X 20

Type	Ref No	Description	Case Quantity
Markers			
	F895062	Scleral Double Ended 3.5mm & 4mm	12 X 20
Needle Holders			
Combination			
	F893122	Foster with R/H Scissors 12.5cm (5")	12 X 20
Micro			
	F895037	Barraquer Curved With Lock 10cm (4")	12 X 20
	F895036	Barraquer Curved Without Lock 10cm (4")	12 X 20
	F895034	Castroviejo Curved With Lock 14cm (5.5")	12 X 20
	F895032	Castroviejo Curved Without Lock 14cm (5.5")	12 X 20
	F895035	Castroviejo Straight With Lock 14cm (5.5")	12 X 20
	F895033	Castroviejo Straight Without Lock 14cm (5.5")	12 X 20
	F895038	Silcock With Lock 14.5cm (5.8")	12 X 20
	F895039	Troutman Curved With Lock 14cm (5.5")	12 X 20
Probes			
Lacrimonal			
	F895083	Double Ended 0000/000 Liebreichs	12 X 20
	F895069	Double Ended 00/0 Liebreichs	12 X 20
	F895040	Double Ended 1/2 Liebreichs	12 X 20
	F895085	Double Ended 3/4 Liebreichs	12 X 20
	F895064	Foster Set of Three 0.6mm, 0.8mm, 1.0mm	12 X 20
Forceps			
Eye			
	F895028	Graefe Iris Toothed 1:2 8.75cm (3.5")	12 X 20
	F893065	Iris Non Toothed 10cm (4")	12 X 20
	F893066	Iris Toothed 10cm (4")	12 X 20
	F895019	Iris Toothed Fine 10cm (4")	12 X 20
	F895015	Moorfields Non Toothed 11.5cm (4.6")	12 X 20
	F895014	St Martins 1:2 Teeth and Platform 0.3mm 9cm (3.6")	12 X 20
	F895024	St Martins Non Toothed 9cm (3.6")	12 X 20

Type	Ref No	Description	Case Quantity
Forceps (Continued)			
	Fixation		
	F895070	Capsulorhexis Forceps Long	12 X 20
	F895016	Curved Non Toothed 11cm (4.4")	12 X 20
	F895029	Jayle 1:2 Teeth 11.25cm (4.5")	12 X 20
	F895017	Serrated Non Toothed Straight 12cm (4.8")	12 X 20
	F895018	Straight 1:2 Teeth 11cm (4.4")	12 X 20
	Tying		
	F895026	Birks (Tennant) Micro Curved 10cm (4")	12 X 20
	F895027	Birks (Tennant) Micro Straight 10cm (4")	12 X 20
	Corneal		
	F895022	St Martins 1:2 Teeth and Platform 0.15mm 9cm (3.6")	12 X 20
	F895014	St Martins 1:2 Teeth and Platform 0.3mm 9cm (3.6")	12 X 20
	F895024	St Martins Non Toothed 9cm (3.6")	12 X 20
	F895092	Colibri 0.3mm Teeth, No Platform, 7.5cm (5")	12 X 20
	Suture		
	F895021	Moorfields 10cm (4")	12 X 20
	F895023	Moorfields 12cm (4.8")	12 X 20
	Jeweller		
	F897001	Angled Smooth No 7 11cm (4.5")	12 X 20
	F895074	Notched Forceps 11cm (4.5")	12 X 20
	F895080	Straight Very Fine Smooth No 5 11cm (4.5")	12 X 20
	F895079	Straight Fine Smooth No 4 11cm (4.5")	12 X 20
	F897002	Straight Medium Smooth No 3 11cm (4.5")	12 X 20
	Diathermy		
	F895025	Bipolar Forceps 5" and Lead 3m	50
	Epilation		
	F893087	Epilation 8.75cm (3.5")	12 X 20
	F893088	Epilation 10cm (4")	12 X 20

Type	Ref No	Description	Case Quantity
Retractors			
	Hand Held		
	F893156	Gillies Skin Hook 18cm (7")	8 X 20
	F895042	Watson Orbital Malleable Set of 4	50
	Lid		
	F895094	Desmarres Eye (Small) 12mm	8 X 10
	F895095	Desmarres Eye (Medium) 14mm	8 X 10
	F895041	Desmarres Eye (Large) 16mm	8 X 10
	F895057	Stallard Small	12 X 20
	F895056	Stallard Large	12 X 20
	F895082	Jaffe (Pair)	12 X 20
	Chalazion		
	F895013	Desmarres 13mm Oval Plate	12 X 20
	F895010	Desmarres 16mm Oval Plate	12 X 20
	F895011	Desmarres 18mm Oval Plate	12 X 20
	F895012	Desmarres 20mm Oval Plate	12 X 20
	F895067	Desmarres 31mm Oval Plate	12 X 20
Scissors			
	F895045	Iris C.O.F 11.25cm (4.5")	12 X 20
	F895046	Iris Straight 11cm (4.5")	12 X 20
	F895047	Iris Straight Micro 11cm (4.5")	12 X 20
	F895049	Strabismus Straight Blunt 11.5cm (4.6")	12 X 20
	F895048	Strabismus Curved Blunt 11.5cm (4.6")	12 X 20
	F896014	Stevens Tenotomy Curved Fine 11cm (4.4")	12 X 20
	F895078	Stevens Tenotomy Straight Fine 11cm (4.4")	12 X 20
	F895052	Westcott Tenotomy Straight Sharp 11cm (4.4")	12 X 20
	F895081	Westcott Tenotomy Straight Blunt 11cm (4.4")	12 X 20
	F895044	Enucleation Curved 11cm (4.4")	12 X 20

Type	Ref No	Description	Case Quantity
------	--------	-------------	---------------

Scissors (Continued)

	Micro		
	F895087	Castroviejo Corneal Micro 7mm Blade 10cm (4")	12 X 20
	F895088	Castroviejo Corneal Micro 9mm Blade 10cm (4")	12 X 20
	F895043	Castroviejo Corneal Micro 12mm Blade 10cm (4")	12 X 20
	F895050	Vannas Corneal 9mm Blade Curved 8cm (3.2")	12 X 20
	F895091	Vannas Corneal 12mm Blade Curved 8cm (3.2")	12 X 20
	F895051	Vannas Corneal 9mm Blade Straight 8cm (3.2")	12 X 20
	F895090	Vannas Corneal 12mm Blade Straight 8cm (3.2")	12 X 20

Sizers

	Orbital		
	F895084	Sizer 14mm	12 X 20
	F895086	Sizer 16mm	12 X 20
	F895053	Sizer 18mm	12 X 20
	F895054	Sizer 20mm	12 X 20
	F895055	Sizer 22mm	12 X 20

Spatulas

	F895089	Elliot	12 X 20
---	---------	--------	---------

Specula

	Self Retaining		
	F895061	Barraquer Open 40 x 26.5mm	12 X 20
	F895058	Barraquer Open 44 x 24mm	12 X 20
	F895075	Kratz Barraquer Open 40 x 33mm Angled	12 X 20
	F895076	Kratz Barraquer Open 40 x 33mm Straight	12 X 20

Ophthalmic Instrument Sets

Sub Tenon's Set

Ref No	F885002	Forceps Dissecting St. Martins Non Toothed 9cm (3.6")	
UOM	20	Tray	
Contains:		Sterile Field	
Scissors Westcott 11cm (4.4")			
Needle Holder Barraquer Micro, Fine, Curved, No Lock 10cm (4")			

Chalazion Set

Ref No	F885003	Scalpel BP Handle With No 11 Blade	
UOM	20	Tray	
Contains:		Sterile Field	
Retractor Desmarres Eyelid 16mm			
Curette D/E 2mm and 3mm			

Enucleation Strabismus Set

Ref No	F885004	Scissors Strabismus Curved Blunt 11.5cm (4.6")	
UOM	20	Elliot's Spatula	
Contains:		Forceps Artery Halstead Mosquito Straight 12.5cm (5")	
Forceps Moorfields Eye 10cm (4")		Scissors Westcott 11cm (4.4")	
Forceps Sponge Holding Rampley 18cm (7.2")		Needle Holder Barraquer Micro, Fine, Curved, No Lock 10cm (4")	
Hook Strabismus 14cm (5.6")		Gallipot 60ml	
Needle Holder Castroviejo Curved with Lock 14cm (5.5")		Tray	
Clip Bulldog Silver 60mm		Sterile Field	
Forceps Dissecting St. Martins Non Toothed 9cm (3.6")			

Ophthalmic Procedure Packs

Ophthalmic Procedure Packs

Phaco Pack

Ref No	F500071
UOM	20
Contains:	
Bowl Polypropylene 250ml	1
Gallipot Polypropylene 60ml	1
Swab Non Woven 10 x 10cm 4Ply	10
Swab Visispears Lint Free	5
Knife Slit 2.75mm Angled Bevel	1
Cannula Hydrodissector Nucleus 27g x 7/8" Sauter	1
Cannula Lacrimal 21g x 1 3/8"	1
Cannula Anterior Chamber Rycroft 30g x 22mm	1
Syringe 3ml Luer Lock	1
Eye Shield	1
Drape Arm Sleeve 55cm	2
Drape Universal Ophthalmic with Full Incise, Pouch and Wick 140 x 140cm	1
Set of Drug Labels	1



Sharps Pad	1
Speculum Barraquer	1
Forceps Dissecting Moorfields 12cm (4.6")	1
Scissors Strabismums Straight Blunt 11.5cm (5.6")	1
Scissors Westcott 11cm (4.4")	1
Forceps Sponge Holding Rampley 18cm (7.2")	1
Forceps Colibri 0.3mm, Teeth, No platform, 7.5 cm (5")	1
Tray	1
Sterile Field	1

Sub Tenon's Pack

Ref No	F500072
UOM	20
Contains:	
Swab Gauze 5 x 5cm 8ply	5
Gallipot 60ml	1
Syringe 5ml Luer Lock	1
Cannula Sub Tenon Curved 19g x 25mm	1
Eye Pad 8 x 6cm Non Woven	1
Needle Hypodermic 21g x 1.5" Green	1
Drape 76 x 76cm with Central Adhesive Fenestration 6 x 8cm	1



Forceps Moorfields Eye 10cm (4")	1
Scissors Westcott 11cm (4.4")	1
Speculum Barraquer	1
Tray	1
Sterile Field	1

Ophthalmic Procedure Packs (Continued)

AMD Pack

Ref No F500073

UOM 20

Contains:

Braunstein Fixed Caliper 3.5 & 4mm	1
Speculum Barraquer Adult	1
Forceps Moorfields Eye 10cm (4")	1
Drape Fenestrated 40 x 40cm with Adhesive Fenestration 6cm	1
Scissors Iris Micro Eye Straight 11.25cm (4.5")	1
Swab Non Woven 7.5 x 7.5cm 4ply	5
Gallipot 60ml	2
Forceps Sponge Holding Rampley 18cm (7.2")	1



Swab Absorbent Eye On A Stick	5
Cotton Tipped Applicator Small Head 15cm Long	3
Marker Scleral D/E 3.5mm & 4mm	1
Tray	1
Sterile Field	1


Oculoplastic Pack

Ref No F500074

UOM 20

Contains:

Bowl Polypropylene 250ml	1
Gallipot Polypropylene 60ml	2
Applicator Cotton Tipped Small Head 15cm Long	10
Swab Gauze XRD 10 x 10cm 16ply	10
Forceps Multipurpose Blue	4
Drape Side Adhesive & Reinforced 90 x 90cm	2
Drape Bottom Adhesive & Reinforced 175 x 180cm	1
Sharps Pad	1
Forceps St Martin Tying 1:2 Teeth & Platform	1
Forceps Sponge Holding Rampley 18cm (7.2")	1
Scissors Westcott 11cm (4.4")	1



Forcep Halstead Mosquito Curved Fine Tips 12.5cm (5")	1
Forceps Desmarres Chalazion 13mm Oval	1
Scalpel BP Handle With No 15 Blade	1
Needle Holder Castroviejo Straight 14cm (5.5")	1
Scissors Stevens Tenotomy Curved 11cm (4.5")	1
Speculum Barraquer	1
Tray	1
Sterile Field	1

numerical full index

Product Code	Page
F894153	89
F894154	89
F894156	92
F894157	92
F894158	81
F894159	81
F894160	82
F894161	83
F894162	82
F894163	82
F894164	82
F894165	82
F894166	82
F894167	82
F895002	99
F895003	99
F895004	99
F895005	99
F895006	81
F895007	99
F895008	99
F895009	99
F895010	99, 102
F895011	99, 102
F895012	99, 102
F895013	99, 102
F895014	100, 101
F895015	100
F895016	101
F895017	101
F895018	101
F895019	100
F895020	99
F895021	101
F895022	101
F895023	101
F895024	100, 101
F895025	101
F895026	101
F895027	101
F895028	100
F895029	101
F895030	99
F895031	99
F895032	100

Product Code	Page
F895033	100
F895034	100
F895035	100
F895036	100
F895037	100
F895038	100
F895039	100
F895040	100
F895041	102
F895042	102
F895043	103
F895044	102
F895045	102
F895046	102
F895047	102
F895048	40, 102
F895049	40, 102
F895050	103
F895051	103
F895052	102
F895053	103
F895054	103
F895055	103
F895056	102
F895057	102
F895058	103
F895061	103
F895062	100
F895064	100
F895066	99
F895067	99, 102
F895068	99
F895069	100
F895070	101
F895071	99
F895072	99
F895073	99
F895074	101
F895075	103
F895076	103
F895077	99
F895078	102
F895079	101
F895080	101
F895081	102

Product Code	Page
F895082	102
F895083	100
F895084	103
F895085	100
F895086	103
F895087	103
F895088	103
F895089	103
F895090	103
F895091	103
F895092	101
F895093	99
F895094	102
F895095	102
F896001	109
F896002	109
F896003	109
F896004	109
F896005	109
F896006	109
F896007	23, 40, 84, 109
F896008	33, 109
F896009	34, 110
F896011	110
F896012	110
F896013	110
F896014	89, 102
F896015	111
F896016	109
F896017	109
F896018	113
F896019	113
F896020	113
F896021	113
F896022	113
F896023	114
F896024	114
F896026	114
F896028	114
F896029	114
F896030	113
F896031	113
F896032	113
F896033	93
F896034	113

numerical full index

Product Code	Page
F896035	113
F896036	113
F896037	113
F896038	113
F896040	113
F896041	113
F896042	113
F896043	113
F896045	109
F896046	124
F896047	109
F896048	109
F896049	109
F896050	109
F896051	109
F896052	110
F896053	110
F896054	111
F896055	111
F896056	109
F896057	109
F896058	111
F896059	109
F897001	101
F897002	101
F897003	111
F897004	123
F898002	68
F898003	68
F898004	69
F898005	69
F898008	69
F898009	69
F898011	68
F898012	68
F898016	69
F898017	68
F898018	68
F898019	68
F898020	68
F898021	68
F898022	68
F898023	68
F898024	68
F898025	68

Product Code	Page
F898026	68
F898027	68
F898028	68
F898029	68
F898030	68
F898031	68
F898032	68
F898033	68
F898034	68
F898035	68
F898036	68
F898037	68
F898038	68
F898039	68
F898040	69
F898041	69
F898042	69
F898043	69
F898044	69
F898045	69
F898046	69
F898047	69
F898048	69
F898100	69
F898101	69

conversion chart

Multiply	By	To Obtain
Millimeters (mm)	0.1	Centimeters (cm)
Millimeters (mm)	0.03937	Inches (in)
Centimeters (cm)	0.3937	Inches (in)
Centimeters (cm)	0.01	Meters (m)
Centimeters (cm)	10	Milimeters (mm)
Inches (in)	2.54	Centimeters (cm)

Millimeters

1 mm 2 mm 3 mm 4 mm 5 mm 6 mm 7 mm 8 mm 9 mm 10 mm

French

2 Fr. 4 Fr. 6 Fr. 8 Fr. 10 Fr. 12 Fr. 14 Fr. 16 Fr. 18 Fr. 20 Fr. 22 Fr.

24 Fr. 26 Fr. 28 Fr. 30 Fr. 32 Fr. 34 Fr. 36 Fr. 38 Fr. 40 Fr.

Inches

1/8" 1/4" 3/8" 1/2" 5/8" 3/4" 7/8" 1"

needle and catheter reference chart

Hypodermic Needle Tubing				French Catheter Sizes			Millimeter/Inch/Centimeter Conversion			
	O.D. In.	I.D. Reg. Wall	I.D. Thin Wall	French Size	O.D. mm	O.D. In.	mm	In.	In.	cm
33	0.008	0.004	0.005	1	0.33	0.013	0.25	0.0098	0.25	0.630
32	0.009	0.004	0.005	2	0.66	0.026	0.50	0.0197	0.50	1.270
31	0.010	0.005	0.006	3	1.00	0.039	0.75	0.0295	0.75	1.900
30	0.012	0.006	0.007	4	1.33	0.053	1	0.0394	1	2.540
29	0.013	0.007	0.008	5	1.66	0.066	2	0.0787	2	5.080
28	0.014	0.007	0.008	6	2.00	0.079	3	0.1181	3	7.620
27	0.016	0.008	0.010	7	2.33	0.092	4	0.1574	4	10.160
26	0.018	0.010	0.012	8	2.66	0.105	5	0.1968	5	12.700
25	0.020	0.010	0.012	9	3.00	0.118	6	0.2362	6	15.240
24	0.022	0.012	0.014	10	3.33	0.131	7	0.2755	7	17.780
23	0.025	0.013	0.015	11	3.67	0.144	8	0.3149	8	20.320
22	0.028	0.016	0.018	12	4.00	0.158	9	0.3543	9	22.860
21	0.032	0.020	0.022	13	4.33	0.170	10	0.3937	10	25.400
20	0.035	0.023	0.025	14	4.67	0.184	11	0.4330	11	27.940
19	0.042	0.027	0.031	15	5.00	0.197	12	0.4724	12	30.480
18	0.050	0.033	0.042	16	5.33	0.210	13	0.5118	13	33.020
17	0.059	0.041	0.046	17	5.67	0.223	14	0.5511	14	35.560
16	0.065	0.047	0.052	18	6.00	0.236	15	0.5905	15	38.100
15	0.072	0.054	0.059	19	6.33	0.249	16	0.6299	16	40.640
14	0.083	0.063	0.071	20	6.67	0.263	17	0.6692	17	43.180
13	0.095	0.071	0.077	21	7.00	0.276	18	0.7086	18	45.720
12	0.109	0.085	0.091	22	7.33	0.288	19	0.7483	19	48.260
11	0.120	0.094	0.100	23	7.67	0.302	20	0.7874	20	50.800
10	0.134	0.106	0.114	24	8.00	0.315	21	0.8267	21	53.340
9	0.148	0.118	0.126	25	8.33	0.328	22	0.8661	22	55.880
8	0.165	0.135	0.143	26	8.67	0.341	23	0.9055	23	58.420
7	0.180	0.150	0.158	27	9.00	0.354	24	0.9448	24	60.960
6	0.203	0.173	0.181	28	9.33	0.367	25	0.9842	25	63.500



 SignatureTM

"Make dealing with Unisurge a standard procedure"



Suite 12, Bunkilla Plaza,
Bracetown Business Park,
Clonee, Dudlin 15,
D15XC64, Republic of Ireland.

Tel: + 353 1 806 8600
Direct: + 353 1 834 4340
Email: info@tpwhelehan.ie
Web: www.whelehansurgical.ie

FOR MORE DETAILS ON OUR COMPREHEN

For Ordering, please contact:

Uniphar Hospital Services,
4045 Kingswood Road,
Citywest Business Park, Dublin 24.
E-Mail: Unipharhospitalservices@uniphar.ie

Tel: + 353 1 468 8456
Fax: + 353 1 466 0782